




EPOXY SINKS	3	DROP-IN	3
		CUP	4
		END AND WALL	4
		ACCESSORIES	5
STAINLESS STEEL SINKS	6	18-GAUGE	6
		20-GAUGE	8
		22-GAUGE	9
		ACCESSORIES	9
POLYMER SINKS	10		
SERVICE FITTINGS	11	BALL VALVES	11
		NEEDLE VALVES	11
		WATER VALVES	12
		LABORATORY FAUCETS	13
		LABORATORY FAUCET	
		ACCESSORIES	14
		STANDARD FAUCETS	15
		STANDARD FAUCET	
		ACCESSORIES	16
ELECTRICAL FITTINGS	17		
APPARATUS ASSEMBLIES	18		
PEGBOARDS	19		
SAFETY STATIONS	20		
METAL SAFETY CABINETS	21		
WOODEN SAFETY CABINETS	23		
SAFETY EQUIPMENT	26	GERMICIDAL CABINET	26
		GOGGLES	26
		SPILL CONTROL	
		STATION	26
		FIRE BLANKET	
		IN CABINET	26





## DROP-IN RECTANGULAR EPOXY SINKS





Model #	Outlet	Dimensions (in)							
		Inside			Outside			Outlet	
		A	B	C	D	E	F	G	H
52D03C	Center	12.0	8.0	5.8	13.6	9.6	6.8	4.0	6.0
52D05	Corner	14.0	10.0	6.2	15.6	11.6	7.2	3.5	3.5
52D10	Corner	16.0	8.0	6.8	17.6	9.6	7.8	4.0	3.5
52D20	Corner	16.0	16.0	7.5	17.6	17.6	8.5	3.5	3.5
52D25	Corner	18.0	15.0	7.9	19.6	16.6	8.9	3.5	3.5
52D30	Corner	18.0	15.0	11.0	19.6	16.6	11.9	3.5	3.5
52A55	Corner	25.0	15.0	4.8	26.6	16.6	5.8	3.5	3.5
52D55	Corner	25.0	15.0	10.0	26.6	16.6	11.0	4.5	4.5
52D65C	Center	35.5	19.5	9.8	37.1	21.2	10.8	9.8	17.8

## DROP-IN HEXAGONAL EPOXY SINK




Model #	Outlet	Inside Dimensions (in)		
		Length	Width	Depth
52DH20C	Center	30.0	26.0	6.9

## DROP-IN ELLIPTICAL EPOXY SINK




Model #	Outlet	Inside Dimensions (in)				
		Rear Length	Front Length	Front Radius	Apex Width	Depth
52D99	Center	23.0	32.1	R26	18.0	10

## ROUND CUP EPOXY SINK



Model #	Dimensions (in)					
	Inside			Outside		
	L	W	D	L	W	D
52TS07	6.0	6.0	6.0	7.4	7.4	11.4


## OVAL CUP EPOXY SINK



Model #	Dimensions (in)					
	Inside			Outside		
	L	W	D	L	W	D
52CS4R	5.7	2.7	3.7	7.5	4.5	8.3
52TS12	8.8	2.7	5.8	10.4	4.5	10


## RECTANGULAR CUP EPOXY SINK

Includes stainless steel screen.



Model #	Dimensions (in)					
	Inside			Outside		
	L	W	D	L	W	D
52TS18	13.6	4.4	5.4	15.2	6.0	9.9


## END AND WALL EPOXY SINK



Model #	Outlet	Dimensions (in)			
		L	Inside W	D	Back Splash
52E46	Corner	22.0	13.6	12.6	3.8


### END STRAINER/ OVERFLOW

With top disc strainer/  
beehive.



**52POE4** 4" Overflow  
**52POE6** 6" Overflow  
**52POE8** 8" Overflow


### OUTLET/ STOPPER



**52SO3R** Outlet  
**52PSS2** Stopper

### UNIVERSAL TRAP


"P" or "S" configuration.



Model #	Inlet	Outlet	Dimensions (in)						
			A	A1	B		B1		C
					MIN	MAX	MIN	MAX	
<b>52P91</b>	1 1/2 NPS	1 1/2 NPT	9 1/8	6	2 1/2	7 1/4	5 3/4	10 1/2	2 1/2

### BOTTLE TRAP


Made of polyolefin  
material. Conforms to  
ASTM #F1412-92.



Model #	Inlet	Outlet	Dimensions (in)		
			A	B	
			MIN	MAX	
<b>52P96</b>	1 1/2 NPS	1 1/2 NPT	3	3	7 1/4

### ANTI-SYPHON BOTTLE TRAP

Anti-syphon action.



Model #	Size	Dimensions (in)			
		A	B		C
			MIN	MAX	
<b>52P95</b>	1 1/2	3	2	4 1/4	7 3/8

## CUSTOM SINGLE BOWL

### 18-GA./304 ALLOY

Specify faucet hole punching required, per illustrations shown on page 8.



Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
52SL1515AGR	15	6	15	12	9	3 1/2
52SL1613AGR	13	7 1/2	16	10	10	3 1/2
52SL1617AGR	17	7 1/2	16	14	10	3 1/2
52SL17519AGR	19	7 1/2	17 1/2	16	11 1/2	3 1/2
52SL1815AGR	15	7 1/2	18	12	12	3 1/2
52SL1921AGR	21	7 1/2	19	18	14	3 1/2

## CUSTOM DEEP SINGLE BOWL

### 18-GA./304 ALLOY

Specify faucet hole punching required, per illustrations shown on page 8.



Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
52SLX1617AGR	17	10	16	14	10	3 1/2
52SLX1815AGR	15	9	18	12	12	3 1/2
52SLX1921AGR	21	10 1/2	19	18	14	3 1/2

## CUSTOM ADA SINGLE BOWL

### 18-GA./304 ALLOY

Specify faucet hole punching required, per illustrations shown on page 8.



Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H specify depth in 1/2" increments (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
52SLADA1613AGR	13	4 1/2 - 6 1/2	16	10	10	3 1/2
52SLADA1617AGR	17	4 1/2 - 6 1/2	16	14	10	3 1/2
52SLADA17519AGR	19	4 1/2 - 6 1/2	17 1/2	16	11 1/2	3 1/2
52SLADA1815AGR	15	4 1/2 - 6 1/2	18	12	12	3 1/2
52SLADA1921AGR	21	4 1/2 - 6 1/2	19	14	18	3 1/2

## CUSTOM DOUBLE BOWL

### 18-GA./304 ALLOY

Specify faucet hole punching required, per illustrations shown on page 8.



Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
52DL1829AGR	29	7 1/2	18	12	12	3 1/2
52DL1925AGR	25	7 1/2	19	10	14	3 1/2
52DL1933AGR	33	7 1/2	19	14	14	3 1/2

## CUSTOM DEEP DOUBLE BOWL

### 18-GA./304 ALLOY

Specify faucet hole punching required, per illustrations shown on page 8.



Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
52DLX1829AGR	29	9	18	12	12	3 1/2
52DLX1933AGR	33	10 1/2	19	14	14	3 1/2

## CUSTOM ADA DOUBLE BOWL

### 18-GA./304 ALLOY

Specify faucet hole punching required, per illustrations shown on page 8.



Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H specify depth in 1/2" increments (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
52DLADA1829AGR	29	4 1/2 - 6 1/2	18	12	12	3 1/2
52DLADA1925AGR	25	4 1/2 - 6 1/2	19	10	14	3 1/2
52DLADA1933AGR	33	4 1/2 - 6 1/2	19	14	14	3 1/2

## CUSTOM SINGLE BOWL

### 20-GA./304 ALLOY

Specify faucet hole punching required, per illustrations shown below.



Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
52SL1515BGR	15	6	15	12	9	3 1/2
52SL1617BGR	17	7	16	14	10	3 1/2
52SL17519BGR	19	7	17 1/2	16	11 1/2	3 1/2
52SL1921BGR	21	7	19	18	14	3 1/2

## CUSTOM DOUBLE BOWL

### 20-GA./304 ALLOY

Specify faucet hole punching required, per illustrations shown below.



Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
52DL1829BGR	29	7	18	12	12	3 1/2
52DL1925BGR	25	7	19	10	14	3 1/2
52DL1933BGR	33	7	19	14	14	3 1/2


### Faucet Hole Punching Options – Single Bowl Sinks

Specify faucet hole punching required by adding the appropriate number code.


- 1 for (1) hole, centered
- 2 for (2) holes on 4" centers
- 3 for (3) holes on 4" centers

For example:  
specify 52SL1515BGR-3 for 52SL1515BGR sink with 3 holes on 4" centers


(1) hole, centered



(2) holes on 4" centers



(3) holes on 4" centers




### Faucet Hole Punching Options – Double Bowl Sinks

Specify faucet hole punching required by adding the appropriate number code.


- 1 for (1) hole, centered
- 2 for (2) holes on 4" centers
- 3 for (3) holes on 4" centers
- 4 for (4) holes on 4" centers

For example:  
specify 52DL1829BGR-3 for 52DL1829BGR sink with 3 holes on 4" centers


(1) hole, centered




(2) holes on 4" centers



(3) holes on 4" centers



(4) holes on 4" centers



## SINGLE BOWL

### 22-GA. STAINLESS STEEL

Specify faucet hole punching required, per illustrations on page 8.



Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
5211516	15	5 1/8	15	12	10	3 1/2
5211719	17	6 1/8	19	14	14	3 1/2

## DOUBLE BOWL

### 22-GA. STAINLESS STEEL

Specify faucet hole punching required, per illustrations on page 8.



Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
5223319	33	6 1/4	19	14	14	3 1/2
5222519	25	6 1/4	19	10	14	3 1/2

## DRAIN FITTING

ADA offset drain. Stainless steel cup strainer with removable basket. 1 1/2" OD chrome plated brass offset.



52JADA35 Drain Fitting

## STRAINER

Stainless steel cup strainer with removable basket. Basket fitted with rubber stopper. 1 1/2" OD chrome plated 4" tailpiece.



52J35 Strainer

## DOUBLE BOWL WASTE CONNECTION


Continuous waste for double compartment sink. 1 1/2" OD chrome plated, with slip fitting at center outlet. Specify dimensions centerline to centerline of drain openings.



52J52S Double Bowl Waste Connection

**RECTANGULAR  
SINGLE BOWL**
**POLYMER  
INTEGRAL  
UNDER-MOUNT**


Available in White, Soft  
White and Creme.



Model #	Bowl Size (in)			Overall Sink Size (in)			Outlet
	A (left to right)	B (front to back)	C (bowl depth)	D (left to right)	E (front to back)	F (depth)	
52K160	20 3/4	15 3/4	9	23	18	9 7/8	3 9/16

**EQUAL  
DOUBLE BOWL**
**POLYMER  
INTEGRAL  
UNDER-MOUNT**


Available in White, Soft  
White and Creme.



Model #	Bowl Size (in)			Overall Sink Size (in)			Outlet
	A (left to right)	B (front to back)	C (bowl depth)	D (left to right)	E (front to back)	F (depth)	
52K130	13 11/16	16 15/16	8 7/8	30 7/8	18 15/16	9 3/4	3 9/16

**OFFSET  
DOUBLE BOWL**
**POLYMER  
INTEGRAL  
UNDER-MOUNT**


Available in White, Soft  
White and Creme.



Model #	Large Bowl Size (in)			Overall Sink Size (in)			Outlets
	A (left to right)	B (front to back)	C (bowl depth)	D (left to right)	E (front to back)	F (depth)	
52K140	17 1/2	16 1/4	9	31 1/2	18 1/4	9 7/8	3 9/16
	Small Bowl Size (in)						
	A <sub>1</sub> (left to right)	B (front to back)	C <sub>1</sub> (bowl depth)				
	10 1/2	16 1/4	6 1/2				

**OVAL  
VANITY**
**POLYMER  
INTEGRAL  
UNDER-MOUNT**

Available in White, Soft  
White and Creme.



Model #	Bowl Size (in)			Overall Sink Size (in)			Outlet
	A (left to right)	B (front to back)	C (bowl depth)	D (left to right)	E (front to back)	F (depth)	
52V130	16 1/2	13	5 3/8	18 1/2	15	7 1/4	1 3/4

## BALL VALVES

- Forged brass with polished chrome finish
- Serrated hose end
- Chrome plated brass ball and molded TFE seals
- On/off control of air, gas and vacuum
- Working pressure up to 75 PSI
- Chrome plated forged lever handle
- 3/8" IPS mounting shank

### BALL VALVE

#### AIR/GAS/VACUUM PANEL MOUNT

Specify index code by using symbol shown below as a suffix to model numbers (e.g., 52200 AIR).



- ┌ 52200 Standard, Single Valve  
 K 52201 Standard, Double Valve  
     90° wye
- ┌ 52203V Vandal Resistant, Single Valve  
 K 52204V Vandal Resistant, Double Valve  
     90° wye

### BALL VALVE

#### AIR/GAS/VACUUM DECK MOUNT

Specify index code by using symbol shown below as a suffix to model numbers (e.g., 52206 AIR).



- 52206 Standard, Single Valve  
 ○ 52207 Standard, Double Valve  
     90°
- 52208 Standard, Double Valve  
     180°
- ⊗ 52209 Standard, Quad Valve
- 52211V Vandal Resistant, Single Valve  
 ○ 52212V Vandal Resistant, Double Valve  
     90°
- 52213V Vandal Resistant, Double Valve  
     180°
- ⊗ 52214V Vandal Resistant, Quad Valve

## NEEDLE VALVES

- Integral serrated hose end
- Floating S.S. needle and integral brass seat
- 3/8" NPT inlet
- Chrome plated forged four-arm handle
- Regulated control for pressures up to 150 PSI

### NEEDLE VALVE

#### AIR/GAS/VACUUM PANEL MOUNT

Specify index code by using symbol shown below as a suffix to model numbers (e.g., 52216 AIR).



- ┌ 52216 Standard, Single Valve  
 K 52217 Standard, Double Valve  
     90° wye

### NEEDLE VALVE

#### AIR/GAS/VACUUM DECK MOUNT

Specify index code by using symbol shown below as a suffix to model numbers (e.g., 52221 AIR).



- 52221 Standard, Single Valve  
 ○ 52222 Standard, Double Valve  
     90°
- 52223 Standard, Double Valve  
     180°
- ⊗ 52224 Standard, Quad Valve

STANDARD INDEXING	INDEX COLOR	LETTER COLOR	SYMBOL
Air	Orange	Black	AIR
Cold Water	Green	White	CW
Deionized Water	White	Black	DI
Distilled Water	White	Black	DW
Gas	Blue	White	GAS
Hot Water	Red	White	HW
Oxygen	Lt. Green	Black	OXY
Vacuum	Yellow	Black	VAC

Specify index code for color-coded discs by using symbols shown as a suffix to model numbers (e.g., 52200 GAS).

**WATER VALVE**

**DECK MOUNT**

Pure water valve. Single outlet straight valve. 3/8" IPS mounting shank. Tin-lined for distilled, deionized or reverse osmosis water.



**52236** Single Valve  
**52237** Double Valve 180°

**WATER VALVE**

**OVERHEAD MOUNT**

Overhead mount water valve assembly. Turret base with inlet and outlet at 90°. 3/8" IPS mounting shank. Serrated hose end. Forged brass handle.



**52238** Single Valve

**WATER VALVE**

**OVERHEAD MOUNT**

Overhead mount double water valve assembly. 3/8" IPS mounting shank. Serrated hose ends. Forged brass handles.



**52239** Double Valve

## LABORATORY FAUCETS


- Faucet bodies: Cast brass with chrome plated finish
- Outlets: 3/8" IPS female with removable ten serration hose end
- Valves: Self-contained renewable with replaceable stainless steel seats
- Faucets provided with brass four-arm handles with color-coded index discs
- Faucets assembled and factory tested

### SINGLE SERVICE FAUCET

#### VACUUM BREAKER DECK MOUNT

Deck mounted for hot or cold water. 6" rigid gooseneck. Integral vacuum breaker. 1/2" IPS with 1/2" OD slipjoints.

Specify index code by using symbol shown on page 11 as a suffix to model numbers (e.g., 52260-6 HW).




<b>52260-6</b>	Standard LH
<b>52261-6</b>	Standard 180°
<b>52262-6</b>	Standard RH
<b>52260-6WB</b>	Standard LH / Wrist Blade
<b>52261-6WB</b>	Standard 180° / Wrist Blade
<b>52262-6WB</b>	Standard RH / Wrist Blade
<b>52264V-6</b>	Vandal Resistant LH
<b>52265V-6</b>	Vandal Resistant 180°
<b>52266V-6</b>	Vandal Resistant RH
<b>52264V-6WB</b>	Vandal Resistant LH / Wrist Blade
<b>52265V-6WB</b>	Vandal Resistant 180° / Wrist Blade
<b>52266V-6WB</b>	Vandal Resistant RH / Wrist Blade

### SINGLE SERVICE FAUCET

#### VACUUM BREAKER PANEL MOUNT

Panel mounted for hot or cold water. 6" rigid/ swing gooseneck. Integral vacuum breaker. 1/2" IPS with 1/2" OD slipjoints.

Specify index code by using symbol shown on page 11 as a suffix to model numbers (e.g., 52268-6 HW).



<b>52268-6</b>	Standard LH
<b>52269-6</b>	Standard RH
<b>52268-6WB</b>	Standard LH / Wrist Blade
<b>52269-6WB</b>	Standard RH / Wrist Blade

## PURE WATER FAUCET

### DECK MOUNT

Deck mounted for distilled, deionized or reverse osmosis water. Body and gooseneck dipped in liquid tin to coat all interior surfaces. Forged brass with polished chrome plated finish. Self-contained renewable unit with replaceable stainless steel seat. Valve has self-closing control. White nylon four-arm handle with color-coded index disc. 6" rigid gooseneck in right hand configuration. 3/8" IPS male shank. Removable serrated hose end.



**52235** Standard

Specify index code by using symbol shown on page 11 as a suffix to model numbers (e.g., 52235 DW).

## MIXING FAUCET

### VACUUM BREAKER DECK MOUNT

Deck mounted for hot and cold water. 6" rigid gooseneck. Integral vacuum breaker. 1 3/16" male shank with two 1/2" OD slipjoints.

Specify index code by using symbol shown on page 11 as a suffix to model numbers (e.g., 52244-6 HW).



<b>52244-6</b>	Standard
<b>52244-6WB</b>	Standard / Wrist Blades
<b>52246V-6</b>	Vandal Resistant
<b>52246V-6WB</b>	Vandal Resistant / Wrist Blades

## MIXING FAUCET

### SWING SPOUT DECK MOUNT

Deck mounted for hot and cold water. 9" swing spout. 1 3/16" male shank with two 1/2" OD slipjoints.

Specify index code by using symbol shown on page 11 as a suffix to model numbers (e.g., 52248-9 HW).



<b>52248-9</b>	Standard
<b>52248-9WB</b>	Standard / Wrist Blades

## DOUBLE SERVICE FAUCET

### DECK MOUNT

Deck mounted for hot or cold water. Two ball valves for either gas, air or vacuum. 6" rigid gooseneck. 3/8" IPS inlets.

Specify index code by using symbol shown on page 11 as a suffix to model numbers (e.g., 52286-6 HW).



**52286-6** Integral Vacuum Breaker

## HEAVY DUTY DOUBLE SERVICE FAUCET

### DECK MOUNT

Deck mounted for hot or cold water. 1/2" heavy wall brass pipe riser. Two ball valves for either gas, air or vacuum. 3/8" IPS inlets.

Specify index code by using symbol shown on page 11 as a suffix to model numbers (e.g., 52291 HW).



**52291** Inline Vacuum Breaker  
**52293V** Vandal Resistant / Inline Vacuum Breaker

## DOUBLE SERVICE FAUCET

### DECK MOUNT

Deck mounted for hot or cold water. One ball valve for either gas, air or vacuum. 6" rigid gooseneck. One plug on opposite side. 3/8" IPS inlets.

Specify index code by using symbol shown on page 11 as a suffix to model numbers (e.g., 52288-6 HW).



**52288-6** Integral Vacuum Breaker

## MIXING FAUCET

### PANEL OR CURB MOUNT

Panel or curb mount for hot and cold water. 8" rigid/swing convertible gooseneck. 3/8" IPS female union inlets. 8" inlet centers.

Specify index code by using symbol shown on page 11 as a suffix to model numbers (e.g., 52277-8 HWCW).



**52277-8** Integral Vacuum Breaker  
**52277-8WB** Integral Vacuum Breaker / Wrist Blades

## 4" WRIST BLADE HANDLE

4" wrist blade handle with index discs for water faucets.

Specify index code by using HW or CW as a suffix to model number (HW for Hot Water; CW for Cold Water (e.g., 522WB HW).



**522WB** 4" Wrist Blade Handle

## VACUUM BREAKER

Faucet vacuum breaker mounts between outlet gooseneck end and removable serrated hose connection.



**522FV** Vacuum Breaker

## AERATOR

Anti-splash aerator can be used in lieu of serrated hose connection.



**522AR** Aerator

**STANDARD FAUCETS**

- Single lever and handle models
- Available models with spray
- Swing spouts and gooseneck designs
- Typical for applications with stainless steel sinks
- Chrome finishes

**SINGLE LEVER  
FAUCET**

Chrome finish. Chrome lever handle. Swivel spout. No spray. 8" centers.



**52311** Single Lever Faucet

**DOUBLE LEVER  
GOOSENECK  
FAUCET**

Chrome finish. Chrome lever handles. No spray. 8" centers.



**52316** Double Lever Gooseneck Faucet

**SINGLE LEVER  
FAUCET WITH  
SPRAY**

Chrome finish. Chrome lever handle. Swivel spout. Integral black side spray. 8" centers.



**52312** Single Lever Faucet with Spray

**DOUBLE LEVER  
GOOSENECK  
FAUCET WITH  
SPRAY**

Chrome finish. Chrome lever handles. White side spray. 8" centers.



**52317** Double Lever Gooseneck Faucet with Spray

**BAR  
FAUCET**

Chrome finish. Chrome lever handles. 4" centers.



**52313** Bar Faucet

**BLADE HANDLE  
VANITY FAUCET**

Chrome finish. Chrome blade handles. No pop-up drain assembly. 4" centers.



**52320** Blade Handle Vanity Faucet

**SINGLE LEVER  
VANITY FAUCET  
WITH POP-UP  
DRAIN**

Chrome finish. Chrome lever handle. With pop-up drain assembly. 4" centers.



**52319** Single Lever Vanity Faucet with Pop-Up Drain

**BLADE HANDLE  
VANITY FAUCET  
WITH POP-UP  
DRAIN**

Chrome finish. Chrome blade handles. With pop-up drain assembly. 4" centers.



**52321** Blade Handle Vanity Faucet with Pop-Up Drain

**DRAIN WITH FLAT  
STRAINER**

Chrome plated brass drain with flat grid strainer. 1 1/2" OD chrome plated tailpiece, 4" long.



**52J15FS** Drain with Flat Strainer

**DRAIN WITH FLAT  
STRAINER**








Chrome plated brass drain with flat strainer. 1 1/2" OD chrome plated brass tailpiece, 4" long.



**52J35FS** Drain with Flat Strainer

## ELECTRICAL FITTINGS

- Pedestal boxes with wrinkle black epoxy finish
- 3/4" IPS mounting shanks
- Grey duplex outlet and grey GFI duplex
- Stainless steel cover plates

		20 AMP, 125V STRAIGHT BLADE (DUPLEX)	20 AMP, 125V STRAIGHT BLADE (GFI)
<b>SINGLE GANG FLUSH MOUNT</b>   (duplex shown)		<b>52366D</b> (1) Duplex 	<b>52366G</b> (1) GFI 
<b>DOUBLE GANG FLUSH MOUNT</b>   (duplexes shown)		<b>52368D</b> (2) Duplexes 	<b>52368G</b> (1) GFI (1) Duplex 
<b>SINGLE GANG SINGLE FACE PEDESTAL</b>  		<b>52370D</b> (1) Duplex   <b>52371D</b> Vandal Resistant (1) Duplex 	<b>52370G</b> (1) GFI   <b>52371G</b> Vandal Resistant (1) GFI 
<b>DOUBLE GANG SINGLE FACE PEDESTAL</b>  		<b>52373D</b> (2) Duplexes   <b>52374D</b> Vandal Resistant (2) Duplexes 	<b>52373G</b> (1) GFI (1) Duplex   <b>52374G</b> Vandal Resistant (1) GFI (1) Duplex 
<b>SINGLE GANG DOUBLE FACE PEDESTAL</b>  		<b>52376D</b> (2) Duplexes   <b>52377D</b> Vandal Resistant (2) Duplexes 	<b>52376G</b> (1) GFI (1) Duplex   <b>52377G</b> Vandal Resistant (1) GFI (1) Duplex 
<b>DOUBLE GANG DOUBLE FACE PEDESTAL</b>  		<b>52379D</b> (4) Duplexes   <b>52380D</b> Vandal Resistant (4) Duplexes 	<b>52379G</b> (1) GFI (3) Duplexes   <b>52380G</b> Vandal Resistant (1) GFI (3) Duplexes 

**3/4" HORIZONTAL  
APPARATUS RODS**

3/4" aluminum horizontal rods with rounded ends.



	WIDTH
<b>522H1-24</b>	24"
<b>522H1-36</b>	36"
<b>522H1-48</b>	48"

**3/4" VERTICAL  
APPARATUS RODS**

3/4" aluminum vertical rods with tapered end.



	HEIGHT
<b>522V1-24</b>	24"
<b>522V1-36</b>	36"
<b>522V1-48</b>	48"

**3/4" ROD  
CLAMP**

Aluminum. For 3/4" diameter horizontal rods.



**522C1** 3/4" Rod Clamp

**3/4" SURFACE  
MOUNT BURETTE  
SOCKET**

Satin aluminum finish rod socket with tapered bore and radiused top for surface mounting on countertop. For 3/4" diameter vertical rods.



**522SB** 3/4" Surface Mount Burette Socket

**3/4" RECESSED  
BURETTE SOCKET**


Satin aluminum finish rod socket with tapered bore and flat top for recess mounting in countertop. For 3/4" diameter vertical rods.



**522RB** 3/4" Recessed Burette Socket

### STAINLESS STEEL PEGBOARD WITH DRIP DEFLECTOR


Includes wall bracket  
and 6" pegs.



MODEL #	SIZE (W x H)	# PEGS
B2430	24" x 30"	32
B3030	30" x 30"	50
B3630	36" x 30"	50

### STAINLESS STEEL PEGBOARD WITH DRIP TROUGH AND CATCH DRAIN

Includes wall bracket  
and 6" pegs.



MODEL #	SIZE (W x H)	# PEGS
V1824	18" x 24"	15
V2424	24" x 24"	20
V2430	24" x 30"	32
V2436	24" x 36"	40

### EPOXY PEGBOARD

1" BLACK



Pegboard with  
Stainless Steel  
Drip Trough

#### PEGBOARD ONLY

MODEL # unfinished back	MODEL # finished back	W	H	# PEGS
52P1824	52P1824FB	18	24	21
52P2424	52P2424FB	24	24	28
52P2430	52P2430FB	24	30	39
52P3030	52P3030FB	30	30	46
52P3630	52P3630FB	36	30	60
52P3636	52P3636FB	36	36	68

#### PEGBOARD WITH STAINLESS STEEL DRIP TROUGH

MODEL # unfinished back	MODEL # finished back	W	H	# PEGS
52P1824T	52P1824TFB	18	24	21
52P2424T	52P2424TFB	24	24	28
52P2430T	52P2430TFB	24	30	39
52P3030T	52P3030TFB	30	30	46
52P3630T	52P3630TFB	36	30	60
52P3636T	52P3636TFB	36	36	68

## EYE/FACE WASH DECK MOUNTED 90° AUTOFLOW™

### APPLICATION:

AutoFlow™ eye/face wash for mounting on counter next to sink. Water flow is activated by swinging the spray head assembly out over sink.

### CONSTRUCTION:

Polished chrome plated brass.

### MOUNTING:

Unit is mounted on right (52FE774) or left (52FE774LH) side of sink.

### SPRAY HEAD ASSEMBLY:

Two FS-Plus™ spray heads. Each head has a "flip top" dust cover, internal flow control and filter to remove impurities from the water flow.

### VALVE:

1/2" IPS brass plug-type valve with O-ring seals. Swinging spray head assembly from storage to operational position opens orifice and activates water flow. Unit remains in operation until spray head assembly is swung back into the storage position.

### STRAINER:

Unit is furnished with in-line strainer to protect valve and spray heads from debris in water line.

### SUPPLY:

1/2" NPT female inlet.


### SIGN:

Furnished with universal identification sign.



**52FE774LH** Eye/Face Wash Deck Mounted 90° AutoFlow™ - Left Mounted

**52FE774** Eye/Face Wash Deck Mounted 90° AutoFlow™ - Right Mounted




52FE774LH

52FE774

## ADA ACCESSIBLE BARRIER FREE EYE/FACE WASH AND SHOWER

- ADA accessible barrier free eye/face wash and shower
- Complies with all requirements for minimum knee clearance, maximum height and reach and minimum distance from obstructions
- 4 spray heads with shallow pattern bowl to enhance accessibility
- Fully assembled and water tested prior to shipment



**52GBF1909**

ADA Accessible Barrier Free Eye/Face Wash and Shower

## FLAMMABLE LIQUID

Before purchasing any flammable liquid storage cabinet, familiarize yourself with all local and state codes from the appropriate agencies having jurisdiction in your area.

### STAND-ALONE/ STACKABLE BASE



- 30-Gallon
- Manual Doors
- Lever Handle

Model #	Description	Exterior Dimensions	Capacity	Color
52121FY	Stand-Alone/Stackable Base Cabinet	43"W x 44"H x 18"D	30 gallons	Yellow

- 30-Gallon
- Self-Closing
- Lever Handle

Model #	Description	Exterior Dimensions	Capacity	Color
52122FY	Stand-Alone/Stackable Base Cabinet	43"W x 44"H x 18"D	30 gallons	Yellow

### STAND-ALONE



- 45-Gallon
- Manual Doors
- Lever Handle

Model #	Description	Exterior Dimensions	Capacity	Color
52125FY	Stand-Alone Cabinet	43"W x 65"H x 18"D	45 gallons	Yellow

- 45-Gallon
- Self-Closing
- Lever Handle

Model #	Description	Exterior Dimensions	Capacity	Color
52126FY	Stand-Alone Cabinet	43"W x 65"H x 18"D	45 gallons	Yellow

### UNDER-COUNTER/ STAND-ALONE



- 22-Gallon
- Manual Doors
- Lever Handle

Model #	Description	Exterior Dimensions	Capacity	Color
52133FY	Under-Counter/Stand-Alone Cabinet	35"W x 35"H x 22"D	22 gallons	Yellow

### UNDER-COUNTER/ STAND-ALONE



- 12-Gallon
- Manual Door
- Lever Handle

Model #	Description	Exterior Dimensions	Capacity	Color
52129FY	Under-Counter/Stand-Alone Cabinet	23"W x 35"H x 18"D	12 gallons	Yellow

## ACID / CORROSIVES

Before purchasing any chemical storage cabinet, familiarize yourself with all local and state codes from the appropriate agencies having jurisdiction in your area.

### STAND-ALONE/ STACKABLE BASE



- 30-Gallon
- Manual Doors
- Lever Handle

Model #	Description	Exterior Dimensions	Capacity	Color
<b>52121AB</b>	Stand-Alone/Stackable Base Cabinet	43"W x 44"H x 18"D	30 gallons	Blue

- 30-Gallon
- Self-Closing
- Lever Handle

Model #	Description	Exterior Dimensions	Capacity	Color
<b>52122AB</b>	Stand-Alone/Stackable Base Cabinet	43"W x 44"H x 18"D	30 gallons	Blue

### STAND-ALONE



- 45-Gallon
- Manual Doors
- Lever Handle

Model #	Description	Exterior Dimensions	Capacity	Color
<b>52125AB</b>	Stand-Alone Cabinet	43"W x 65"H x 18"D	45 gallons	Blue

- 45-Gallon
- Self-Closing
- Lever Handle

Model #	Description	Exterior Dimensions	Capacity	Color
<b>52126AB</b>	Stand-Alone Cabinet	43"W x 65"H x 18"D	45 gallons	Blue

### UNDER-COUNTER/ STAND-ALONE



- 22-Gallon
- Manual Doors
- Lever Handle

Model #	Description	Exterior Dimensions	Capacity	Color
<b>52133AB</b>	Under-Counter/Stand-Alone Cabinet	35"W x 35"H x 22"D	22 gallons	Blue

## FLAMMABLES / ACID

Before purchasing any flammable liquid or chemical storage cabinet, familiarize yourself with all local and state codes from the appropriate agencies having jurisdiction in your area.

### BENCH FLAMMABLES CABINET



- Constructed of 1" multi-ply, exterior-grade plywood.
- Meets all applicable OSHA and NFPA #30 standards for the storage of flammable liquids.
- Convenient and safe storage of small quantities on the laboratory bench.
- Cabinet floor constructed as a 2" liquid-tight trough.
- One fully-adjustable shelf.

Model #	Description	Exterior Dimensions	Capacity	Color
SC8021	Bench Flammables Cabinet	16¾"W x 21¼"H x 15¾"D	4 gallons	Yellow

### BENCH ACID CABINET



- Wood construction—no metal—no rust or corrosion.
- Constructed of 1" multi-ply, exterior-grade plywood.
- Sized for convenient storage of corrosive chemicals on the laboratory bench.
- Painted with multiple coats of chemical-resistant paint.
- Cabinet floor constructed as a 2" liquid-tight trough.


Model #	Description	Exterior Dimensions	Capacity	Color
SC8063	Bench Corrosives Cabinet	16¾"W x 21¼"H x 15¾"D	5 x 2.5 liter 16 x 500 mL	Blue



### ACID CABINET CONSTRUCTION



Unique, entirely wooden hinge assembly eliminates common problem of hinge deterioration.



Unique polypropylene lock assembly. Padlock and two keys furnished with cabinet.



Top recessed to include a liquid-tight polypropylene tray. Tray will hold 11 liters (3 gallons) of spilled liquid. Use this top tray as a dispensing area.



Cabinet floor constructed as a 2" liquid-tight trough to contain spillage.

## FLAMMABLES / ACID

Before purchasing any flammable liquid or chemical storage cabinet, familiarize yourself with all local and state codes from the appropriate agencies having jurisdiction in your area.

### STAK-A-CAB SAFETY STORAGE CABINETS



SC1460, SC2060, SC1860 Stak-A-Cabs  
with SC1861 Stand

- Constructed of 1" thick, multi-ply plywood and finished with chemical-resistant paint.
- Organize all laboratory chemicals in a space no larger than a standard flammables cabinet.
- Customize a chemical storage center for special storage needs.
- All cabinets constructed to meet and exceed all applicable NFPA and OSHA standards.

### STAK-A-CAB CABINETS

All units 47"W x 18"H x 18"D

Model #	Description	Capacity	Color
SC1360	Stak-A-Cab Acid Cabinet	22 x 2.5 Liter	Blue
SC1460	Stak-A-Cab Corrosives Cabinet	22 x 2.5 Liter	Blue
SC1860	Stak-A-Cab Flammables Cabinet	16 x 1 Gallon	Yellow
SC2060	Stak-A-Cab Combination Acid/Flammables Cabinet	Flammable: 8 x 1 Gallon Corrosive: 10 x 2.5 Liter	Gray

### STAK-A-CAB STANDS

All units 46½"W x 45⅞"H x 17¾"D

Model #	Description	Color
SC1461	Stak-A-Cab Floor Stand for Acid Cabinet	Blue
SC1861	Stak-A-Cab Floor Stand for Flammables Cabinet	Yellow
SC2061	Stak-A-Cab Floor Stand for Acid/Flammables Cabinet	Gray

### JUMBO STACKING CABINETS



- Constructed of 1" thick, multi-ply, exterior-grade plywood.
- Finished with easily identified chemical-resistant paint.
- Organize both flammable and corrosive liquids in a space no larger than a standard cabinet.
- Individual cabinets may be placed under the work surface.
- Design as a stand-alone or stack up to two units high.
- All cabinets meet or exceed all applicable OSHA and NFPA standards.
- All cabinets have a liquid-tight, 2" trough to contain accidental spills.

All units 30"W x 32⅝"H x 18½"D

Model #	Description	Capacity	Color
SC8080	Jumbo Stacking Flammables Cabinet	18 x 1 Gallon	Yellow
SC8081	Jumbo Stacking Acid Cabinet	24 x 2.5 Liter	Blue
SC8082	Jumbo Stacking Corrosives Cabinet	24 x 2.5 Liter	Blue

## FLAMMABLES / ACID

Before purchasing any flammable liquid or chemical storage cabinet, familiarize yourself with all local and state codes from the appropriate agencies having jurisdiction in your area.

### COMBINATION ACID/ FLAMMABLES CABINETS



#### ACID SIDE:

- Wooden hinge assembly.
- Polypropylene hasp assembly.
- 10" corrosive label for easy identification.
- Cabinet floor constructed as 2" liquid-tight trough.
- Color: Light Gray

#### FLAMMABLES SIDE:

- Continuous metal hinge.
- Metal locking assembly.
- Adjustable shelf.
- 10" flammable label for easy identification.
- Cabinet floor constructed as 2" liquid-tight trough.
- Color: Light Gray

Model #	Description	Exterior Dimensions	Capacity	
			Corrosive	Flammable
SC2036	Combination Acid/Flammables Cabinet, Partially Lined	35"W x 35½"H x 22"D	18 x 2.5 Liter	12 x 1 Gallon
SC2048	Combination Acid/Flammables Cabinet, Partially Lined	47"W x 35½"H x 22"D	24 x 2.5 Liter	16 x 1 Gallon

### ACID CABINETS



- Constructed of 1" thick, high density, multi-ply plywood.
- Finished with multiple coats of chemical-resistant paint.
- Cabinets have a liquid-tight 2" trough to contain accidental spills.
- Fasteners used are plastic—even the few screws used are nylon.
- Wooden door handle.
- Word "ACID" printed in 6" red letters.
- 4 leveling feet included.

Model #	Description	Exterior Dimensions	Bottle Capacity
SC8051	Acid Cabinet, Interior Fully Lined with Polypropylene (floor, interior walls, interior top, and top tray)	31"W x 36⅝"H x 20"D	30 x 2.5 Liter or 90 x 500 mL
SC8041	Acid Cabinet, Interior Partially Lined with Polypropylene (floor, shelf, and top tray)	31"W x 36⅝"H x 20"D	30 x 2.5 Liter or 90 x 500 mL

### FLAMMABLES CABINETS



Model #	Description	Approx. Capacity	Interior Dimensions	Exterior Dimensions
SC7021	Floor Flammables Cabinet, 1 Adjustable Shelf	18 Gallon	217⁄8"W x 34"H x 21"D	237⁄8"W x 36⅝"H x 237⁄8"D
SC7131	45-Gallon Flammables Cabinet, 2 Adjustable Shelves	45 Gallon	41"W x 42"H x 151⁄8"D	43"W x 44⅝"H x 18"D

**GERMICIDAL  
CABINET**



For sanitizing and storage of safety glasses or goggles. Holds 48 glasses or 40 goggles. Cycle controlled by timer. Pilot light confirms UV lamp in use. Shuts off automatically if double doors are open. Doors interlocked with tamper-resistant latches. Wall- or shelf-mountable with 7' grounded electrical cord. 24-gauge white enameled lockable steel cabinet. Germicidal lamp. No assembly required. Goggles/Glasses not included. Dimensions: 24 1/2"W x 32"H x 9 1/2"D.

**52192** Germicidal Cabinet

**ADVANTAGE  
GOOGLE**



UV-absorbing polycarbonate lens. Tinted flexible body. Adjustable headband. Indirect vent. Fits over most prescription eyewear. Replaceable lens. Meets ANSI Z87.1-2003 Standards. CUL Certified to meet CAN/CSA-Z94.3 Standards. (direct vent model depicted at left)

**52194** Advantage Goggle

**CHEMICAL SPLASH  
GOOGLE**



Clear UV-absorbing polycarbonate lens. Soft, clear goggle body. Indirect vents top and bottom. Fits over most prescription eyewear. Meets ANSI Z87.1-2003 Standards. CUL Certified to meet CAN/CSA-Z94.3 Standards.

**52196** Chemical Splash Goggle

**SPILL CONTROL  
STATION**



For easy clean-up of emergency spills. Universal spill control pillows absorb and contain all types of acid, caustic, solvent or oil. Special pillows for hydrofluoric acid. Includes:

- 10 standard spill pillows
- 15 sorbent pads
- 2 pairs of gloves
- 15 hazardous waste bags and ties
- 2 hydrofluoric spill pillows
- 2 Total Body Saranex "Tyvek" coveralls
- 2 chemical splash goggles

**52186** Spill Control Station

**FIRE BLANKET  
IN CABINET**



Fire-retardant 65"x84" wool blanket. Useful for general first aid and other emergencies. Heavy 22-gauge red steel cabinet with drop down door. Colorful "Stop, Drop and Roll" technique graphics.

**52188** Fire Blanket in Cabinet