

COUNTERTOPS

AP04 SERIES: PHENOLIC	3	BLACK ACID RESISTANT FACE
AP06 SERIES: PHENOLIC	4	DECORATIVE ACID RESISTANT - COLOR THROUGHOUT
ES00 SERIES: EPOXY	5	SOLID EPOXY RESIN
ES02 SERIES: EPOXY	6	SOLID EPOXY RESIN - MARINE EDGE
LS05 SERIES: LAMINATE	7	HIGH PRESSURE LAMINATE - 3MM EDGE
LS07 SERIES: LAMINATE	8	HIGH PRESSURE LAMINATE - SELF EDGE
LS25 SERIES: ACID RESISTANT LAMINATE	9	HIGH PRESSURE ACID RESISTANT LAMINATE - 3MM EDGE
LS26 SERIES: ACID RESISTANT LAMINATE	10	HIGH PRESSURE ACID RESISTANT LAMINATE - SELF EDGE
SS00 SERIES: STAINLESS STEEL	11	STAINLESS STEEL - FORMED EDGE DESIGN
SS01 SERIES: STAINLESS STEEL	12	STAINLESS STEEL - MARINE EDGE DESIGN
WD00 SERIES: HARDWOOD	13	HARD MAPLE LUMBER

WINDOWSILLS

AP10 SERIES: PHENOLIC	3	BLACK ACID RESISTANT FACE
AP12 SERIES: PHENOLIC	4	DECORATIVE ACID RESISTANT - COLOR THROUGHOUT
ES10 SERIES: EPOXY	5	SOLID EPOXY RESIN

REV. 5/7/21

©2021 STEVENS INDUSTRIES, INC.
ALL RIGHTS RESERVED.
Due to continuing product
development, Stevens reserves
the right to change specifications
without notice.

AP04 SERIES: PHENOLIC

- Black solid phenolic core with black acid resistant face
- Bevel edges, polished for consistent edge look
- Applied backsplash and endsplash
- Drip groove

Available Color:

Black


AT = Phenolic Top
AB = Phenolic Backsplash
AE = Phenolic Endsplash

Series/Core Code:


AP04S Solid

AP04S available countertop thicknesses:
19mm (3/4")
25mm (1")

Shown: AP04S1919



Shown: AP04S2525



AP10 SERIES WINDOWSILLS

AP10 SERIES: PHENOLIC

- Black solid phenolic core with black acid resistant face
- Bevel edge, polished for consistent edge look

Available Color:


Black

AT = Phenolic Top


Series/Core Code:

AP10S Solid

AP10S available windowsill thicknesses:
19mm (3/4")
25mm (1")



Shown: AP10S19



Shown: AP10S25

AP06 SERIES: PHENOLIC

- Decorative acid resistant phenolic, color throughout
- Bevel edges, polished for consistent edge look
- Applied backsplash and endsplash
- Drip groove

Available Colors:

Grey, Sandstone

*Other colors may be available upon request.
Special pricing and quantities may apply.*


AT = Phenolic Top
AB = Phenolic Backsplash
AE = Phenolic Endsplash

Series/Core Code:


AP06S Solid

AP06S available countertop thicknesses:
19mm (3/4")
25mm (1")

Shown: AP06S1919



Shown: AP06S2525



API2 SERIES WINDOWSILLS

API2 SERIES: PHENOLIC

- Decorative acid resistant phenolic, color throughout
- Bevel edge, polished for consistent edge look

Available Colors:

Grey, Sandstone

*Other colors may be available upon request.
Special pricing and quantities may apply.*


AT = Phenolic Top

Series/Core Code:


API2S Solid

API2S available windowsill thicknesses:
19mm (3/4")
25mm (1")

Shown: API2S19



Shown: API2S25



ES00 SERIES: EPOXY

- Solid epoxy resin
- Bevel edges with smooth edge finish
- Applied backsplash and endsplash
- Drip groove
- Butt joints and attachment with field applied epoxy sealant and adhesives
- Epoxy resin coordinated sinks where specified
- Factory fabricated cutouts and drillings

ET = Epoxy Top
EB = Epoxy Backsplash
EE = Epoxy Endsplash


Series/Color Codes:

ES00BL Black
ES00GR Grey
ES00SA Sandstone


*Other colors may be available upon request.
Special pricing and quantities may apply.*

ES00 available countertop thicknesses:
19mm (3/4")
25mm (1")

Shown: ES00BL1919



Shown: ES00BL2525



ES10 SERIES WINDOWSILLS

ES10 SERIES: EPOXY

- Solid epoxy resin
- Bevel edge with smooth edge finish


ET = Epoxy Top

Series/Color Codes:


ES10BL Black
ES10GR Grey
ES10SA Sandstone

*Other colors may be available upon request.
Special pricing and quantities may apply.*

ES10 available windowsill thicknesses:
19mm (3/4")
25mm (1")



Shown: ES10BL19



Shown: ES10BL25

ES02 SERIES: EPOXY

- Solid epoxy resin
- Lipped marine edge
- Bevel edges with smooth edge finish
- Applied backsplash and endsplash
- Butt joints and attachment with field applied epoxy sealant and adhesives
- Epoxy resin coordinated sinks where specified
- Factory fabricated cutouts and drillings

ET = Epoxy Top
EB = Epoxy Backsplash
EE = Epoxy Endsplash


Series/Color Codes:

ES02BL Black
ES02GR Grey
ES02SA Sandstone

*Other colors may be available upon request.
Special pricing and quantities may apply.*

ES02 available countertop thickness:
25mm (1")

Shown: ES02BL2519



LS05 SERIES: LAMINATE

- High pressure laminate
- 3mm edge top
- Applied 3mm edge backsplash and endsplash

Available Colors:

Black, Grey, Sandstone

Other colors may be available upon request.

Special pricing and quantities may apply.

LT = Laminate Top
LB = Laminate Backsplash
LE = Laminate Endsplash


Series/Core Codes:

- LS05C TSCA particleboard
- LS05M Moisture resistant/NAF particleboard
- LS05P Plywood


LS05C, LS05M available countertop thicknesses:
19mm (3/4")
25mm (1")

LS05P available countertop thickness:
25mm (1")

Shown:
LS05C1919



Shown:
LS05C2525



LS07 SERIES: LAMINATE

- High pressure laminate
- Self edge top
- Applied self edge backsplash and endsplash

Available Colors:

Black, Grey, Sandstone

Other colors may be available upon request.

Special pricing and quantities may apply.

LT = Laminate Top
LB = Laminate Backsplash
LE = Laminate Endsplash


Series/Core Codes:

- LS07C TSCA particleboard
LS07M Moisture resistant/NAF particleboard
LS07P Plywood


LS07C, LS07M available countertop thicknesses:
19mm (3/4")
25mm (1")

LS07P available countertop thickness:
25mm (1")

**Shown:
LS07C1919**



**Shown:
LS07C2525**



LS25 SERIES: ACID RESISTANT LAMINATE

- Acid resistant face high pressure laminate
- 3mm edge top
- Applied 3mm edge backsplash and endsplash

Available Colors:

Black, Grey, Sandstone

Other colors may be available upon request.

Special pricing and quantities may apply.

LT = Laminate Top
LB = Laminate Backsplash
LE = Laminate Endsplash


Series/Core Codes:

- LS25C TSCA particleboard
- LS25M Moisture resistant/NAF particleboard
- LS25P Plywood


LS25C, LS25M available countertop thicknesses:
19mm (3/4")
25mm (1")

LS25P available countertop thickness:
25mm (1")

Shown:
LS25C1919



Shown:
LS25C2525



LS26 SERIES: ACID RESISTANT LAMINATE

- Acid resistant face high pressure laminate
- Self edge top
- Applied self edge backsplash and endsplash

Available Colors:

Black, Grey, Sandstone

Other colors may be available upon request.

Special pricing and quantities may apply.

LT = Laminate Top
LB = Laminate Backsplash
LE = Laminate Endsplash


Series/Core Codes:

- LS26C TSCA particleboard
- LS26M Moisture resistant/NAF particleboard
- LS26P Plywood


LS26C, LS26M available countertop thicknesses:
19mm (3/4")
25mm (1")

LS26P available countertop thickness:
25mm (1")

Shown:
LS26C1919



Shown:
LS26C2525



SS00 SERIES: STAINLESS STEEL

- Formed edge design with integral cove
- Fabricated from 16 ga. 304 stainless steel
- Formed edges and bends with welded and polished corners
- Can be fabricated with or without backsplash or endsplash


Series/Core Codes:

SS00R Tube Rail

SS00R available countertop thickness:
29mm (1 1/8")

ST = Stainless Steel Top

Shown:
SS00R2925



Standard Welded-in Integral 16 Ga. Stainless Steel Sink Bowl Sizes*

Left to Right	Front to Back	Depth
10"	14"	8"
10"	14"	10"
16"	18"	12"
16"	20"	12"
18"	18"	12"
18"	24"	12"
20"	20"	12"
24"	24"	12"

*Custom bowl fabrication also available

SS01 SERIES: STAINLESS STEEL

- Marine edge design with integral cove
- Fabricated from 16 ga. 304 stainless steel
- Formed edges and bends with welded and polished corners
- Can be fabricated with or without backsplash or endsplash


ST = Stainless Steel Top

Series/Core Codes:

SS01R Tube Rail

SS01R available countertop thickness:
38mm (1 1/2")

Shown:
SS01R3825



Standard Welded-in Integral 16 Ga. Stainless Steel Sink Bowl Sizes*

Left to Right	Front to Back	Depth
10"	14"	8"
10"	14"	10"
16"	18"	12"
16"	20"	12"
18"	18"	12"
18"	24"	12"
20"	20"	12"
24"	24"	12"

*Custom bowl fabrication also available

WD00 SERIES: HARDWOOD

- Fabricated from laminated kiln dried hard maple lumber
- Industrial grade with butt joints and edge grain fabrication (butcher block) design
- Edges sanded smooth with top and bottom corners broken
- Faces and edges sealed with varnish finish

WT = Hardwood Top
WB = Hardwood Backsplash
WE = Hardwood Endsplash

Series/Core Code:

WD00S Solid hard maple

WD00S available countertop thickness:
38mm (1 1/2")

Shown: WD00S3819

