

84H

3

SHELF	3
SHELF/DRAWER	4
DISPLAY	5
DISPLAY/DRAWER	7
EXPANDED METAL	8
CUBICLE	9
CUBICLE/SHELF	10
CORNER	11
TRAY STORAGE	11
BIN STORAGE	14
FILE DRAWER	15
SMOCK STORAGE	16
WARDROBE STORAGE	17
CHEMICAL STORAGE	19
APPARATUS ROD STORAGE	19
TORSO/SKELETON STORAGE	20
MICROSCOPE STORAGE	21

90H

22

SHELF	22
SHELF/DRAWER	23
DISPLAY	24
DISPLAY/DRAWER	26
EXPANDED METAL	27
CUBICLE	28
CUBICLE/SHELF	29
CORNER	30
SMOCK STORAGE	30
WARDROBE STORAGE	31
APPARATUS ROD STORAGE	33
SAFETY STATION	33
MICROSCOPE STORAGE	34

STACK-ON

35

DESIGN DETAILS	35
17H	36
23H	37
29H	39
35H	41

INSTALLATION ACCESSORIES

43

SOFFIT TRIM	43
VERTICAL SCRIBE	43
FILLERS	43
CLOSURES	44
SOFFIT FASCIAS	45
SOFFIT EXTENSION	45
TOPS	46
OVERLAY BACKS	46
FILE DRAWERS	47

REV. 8/10/21

©2021 STEVENS INDUSTRIES, INC.
ALL RIGHTS RESERVED.
Due to continuing product
development, Stevens reserves
the right to change specifications
without notice.



**ADVANTAGE
SCIENTIFIC**
BY STEVENS

TALL

SHELF


47101

W : 12, 15, 18, 21, 24, 30, 36
H : 84
D : 14, 17, 19, 23

- Fixed Shelf
- Adjustable Shelves: 4

SHELF


47129

W : 30, 36
H : 84
D : 14, 17, 19, 23

- Doors
- Fixed Shelf
- Adjustable Shelves: 4


SHELF


47120

W : 12, 15, 18, 21, 24
H : 84
D : 14, 17, 19, 23

- Right Hinge Door
- Fixed Shelf
- Adjustable Shelves: 4

SHELF


47126

W : 30, 36
H : 84
D : 14, 17, 19, 23

- Upper/Lower Doors
- Adjustable Shelves: 4

SHELF


47121

W : 12, 15, 18, 21, 24
H : 84
D : 14, 17, 19, 23

- Left Hinge Door
- Fixed Shelf
- Adjustable Shelves: 4

SHELF


47138

W : 42, 48
H : 84
D : 14, 17, 19, 23

- Center Divider
- Adjustable Shelves: 10

SHELF


47181

W : 42, 48
H : 84
D : 14, 17, 19, 23

- Doors
- Center Divider
- Adjustable Shelves: 10

SHELF/
DRAWER

47335

W : 30, 36, 42, 48
H : 84
D : 17, 19, 23

- Doors
- Center Divider
- Adjustable Shelves: 6
- Drawers: 3

SHELF


47149

W : 42, 48
H : 84
D : 17, 19, 23

- Sliding Doors
- Center Divider
- Adjustable Shelves: 10

SHELF/
DRAWER

47374

W : 30, 36, 42, 48
H : 84
D : 17, 19, 23

- Doors
- Center Divider
- Adjustable Shelves: 6
- Drawers: 4

SHELF


47150

W : 42, 48
H : 84
D : 17, 19, 23

- Upper/Lower Sliding Doors
- Center Dividers
- Adjustable Shelves: 8

DISPLAY

47130

W : 15, 18, 21, 24
H : 84
D : 17, 19, 23

- Right Hinge Display Door
- Fixed Shelf
- Adjustable Shelves: 4

DISPLAY

47132

W : 15, 18, 21, 24
H : 84
D : 17, 19, 23

- Right Hinge Display Door
- Fixed Shelf
- Adjustable Shelves: 4

DISPLAY

47131

W : 15, 18, 21, 24
H : 84
D : 17, 19, 23

- Left Hinge Display Door
- Fixed Shelf
- Adjustable Shelves: 4

DISPLAY

47133

W : 15, 18, 21, 24
H : 84
D : 17, 19, 23


- Left Hinge Display Door
- Fixed Shelf
- Adjustable Shelves: 4

DISPLAY

47137

W : 30, 36, 42, 48
H : 84
D : 17, 19, 23

- Display Doors
- Center Divider
- Adjustable Shelves: 10

DISPLAY

47140

W : 30, 36, 42, 48
H : 84
D : 17, 19, 23

- Display Doors
- Center Divider
- Adjustable Shelves: 10



DISPLAY



47182

W : 15, 18, 21, 24
H : 84
D : 17, 19, 23

- Right Hinge Upper/Lower Display Doors
- Adjustable Shelves: 4

DISPLAY



47135

W : 30, 36
H : 84
D : 17, 19, 23

- Upper Display Doors
- Lower Doors
- Adjustable Shelves: 4

DISPLAY



47183

W : 15, 18, 21, 24
H : 84
D : 17, 19, 23

- Left Hinge Upper/Lower Display Doors
- Adjustable Shelves: 4

DISPLAY



47154

W : 30, 36, 42, 48
H : 84
D : 17, 19, 23

- Upper Sliding Display Doors
- Lower Sliding Doors
- Center Dividers
- Adjustable Shelves: 8

DISPLAY



47184

W : 30, 36
H : 84
D : 17, 19, 23

- Upper/Lower Display Doors
- Adjustable Shelves: 4

DISPLAY



47158

W : 30, 36, 42, 48
H : 84
D : 17, 19, 23

- Upper Sliding Glass Doors
- Lower Doors
- Center Dividers
- Adjustable Shelves: 8

**DISPLAY/
DRAWER**

47339

W : 30, 36, 42, 48
H : 84
D : 17, 19, 23

- Display Doors
- Center Divider
- Adjustable Shelves: 6
- Drawers: 3

**DISPLAY/
DRAWER**

47379

W : 30, 36, 42, 48
H : 84
D : 17, 19, 23

- Display Doors
- Center Divider
- Adjustable Shelves: 6
- Drawers: 4

**DISPLAY/
DRAWER**

47342

W : 30, 36, 42, 48
H : 84
D : 17, 19, 23

- Sliding Display Doors
- Center Divider
- Adjustable Shelves: 6
- Drawers: 3

**DISPLAY/
DRAWER**

47378

W : 30, 36, 42, 48
H : 84
D : 17, 19, 23

- Sliding Display Doors
- Center Divider
- Adjustable Shelves: 6
- Drawers: 4

**DISPLAY/
DRAWER**

47343

W : 30, 36, 42, 48
H : 84
D : 17, 19, 23

- Sliding Glass Doors
- Center Divider
- Adjustable Shelves: 6
- Drawers: 3

**DISPLAY/
DRAWER**

47383

W : 30, 36, 42, 48
H : 84
D : 17, 19, 23

- Sliding Glass Doors
- Center Divider
- Adjustable Shelves: 6
- Drawers: 4



**EXPANDED
METAL**




47106

W : 15, 18, 21, 24
H : 84
D : 17, 19, 23

- Right Hinge Framed Expanded Metal Door
- Adjustable Shelves: 5
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**



47109

W : 30, 36, 42, 48
H : 84
D : 17, 19, 23

- Framed Expanded Metal Doors
- Center Divider
- Adjustable Shelves: 10
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**




47107

W : 15, 18, 21, 24
H : 84
D : 17, 19, 23

- Left Hinge Framed Expanded Metal Door
- Adjustable Shelves: 5
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**



47108

W : 30, 36
H : 84
D : 17, 19, 23

- Framed Expanded Metal Doors
- Adjustable Shelves: 5
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

CUBICLE


47213

W : 21, 24
H : 84
D : 17, 19, 23

- 3/4" Adjustable Cubicle Shelves: 12
- Compartments: 14
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE


47216

W : 30, 36, 42, 48
H : 84
D : 17, 19, 23

- 3/4" Adjustable Cubicle Shelves: 18
- Compartments: 21
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W), 13" (42W), 15" (48W)

CUBICLE


47214

W : 21, 24
H : 84
D : 17, 19, 23

- Right Hinge Door
- 3/4" Adjustable Cubicle Shelves: 12
- Compartments: 14
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE


47217

W : 30, 36, 42, 48
H : 84
D : 17, 19, 23

- Doors
- 3/4" Adjustable Cubicle Shelves: 18
- Compartments: 21
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W), 13" (42W), 15" (48W)

CUBICLE


47215

W : 21, 24
H : 84
D : 17, 19, 23


- Left Hinge Door
- 3/4" Adjustable Cubicle Shelves: 12
- Compartments: 14
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE


47218


W : 36, 42, 48
H : 84
D : 17, 19, 23

- 3/4" Adjustable Cubicle Shelves: 24
- Compartments: 28
- Compartment Spacing: 3"H min.
- Compartment Widths: 8" (36W), 9 1/2" (42W), 11" (48W)

CUBICLE

47219


W : 36, 42, 48
H : 84
D : 17, 19, 23

- Doors
- 3/4" Adjustable Cubicle Shelves: 24
- Compartments: 28
- Compartment Spacing: 3"H min.
- Compartment Widths: 8" (36W), 9 1/2" (42W), 11" (48W)

**CUBICLE/
SHELF**

47247


W : 36
H : 84
D : 17, 19, 23

- Doors
- Adjustable Shelves: 3
- 3/4" Adjustable Cubicle Shelves: 12
- Compartments: 16
- Compartment Spacing: 3"H min.
- Compartment Width: 8"

**CUBICLE/
SHELF**

47232

W : 30, 36
H : 84
D : 17, 19, 23

- Doors
- Adjustable Shelves: 3
- 3/4" Adjustable Cubicle Shelves: 9
- Compartments: 12
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W)

**CUBICLE/
SHELF**

47246

W : 36
H : 84
D : 17, 19, 23

- Upper/Lower Doors
- Adjustable Shelves: 3
- 3/4" Adjustable Cubicle Shelves: 12
- Compartments: 16
- Compartment Spacing: 3"H min.
- Compartment Width: 8"

**CUBICLE/
SHELF**

47231

W : 30, 36
H : 84
D : 17, 19, 23


- Upper/Lower Doors
- Adjustable Shelves: 3
- 3/4" Adjustable Cubicle Shelves: 9
- Compartments: 12
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W)



SA902: 10W x 3 1/2H x 19D

84H
TALL

**BLIND
CORNER**



47200

W : 33, 36, 39
H : 84
D : 17

W : 36, 39, 42
H : 84
D : 19

W : 39, 42, 45, 48
H : 84
D : 23

- Right Hinge Door
- Fixed Shelf
- Adjustable Shelves: 4

**TRAY
STORAGE**




47263

W : 24
H : 84
D : 23

- SA902 Molded Trays: 32

**BLIND
CORNER**



47201


W : 33, 36, 39
H : 84
D : 17

W : 36, 39, 42
H : 84
D : 19

W : 39, 42, 45, 48
H : 84
D : 23

- Left Hinge Door
- Fixed Shelf
- Adjustable Shelves: 4

**TRAY
STORAGE**




47264

W : 24
H : 84
D : 23

- Right Hinge Door
- SA902 Molded Trays: 32

**TRAY
STORAGE**



47265

W : 24
H : 84
D : 23

- Left Hinge Door
- SA902 Molded Trays: 32




SA901: 14 1/4W x 3 1/2H x 19D



SA902: 10W x 3 1/2H x 19D

84H
TALL

**TRAY
STORAGE**



47270

W : 36
H : 84
D : 23

- SA902 Molded Trays: 48

**TRAY
STORAGE**




47274

W : 48
H : 84
D : 23

- SA901 Molded Trays: 48

**TRAY
STORAGE**




47271

W : 36
H : 84
D : 23

- Doors
- SA902 Molded Trays: 48

**TRAY
STORAGE**




47275

W : 48
H : 84
D : 23

- Doors
- SA901 Molded Trays: 48

**TRAY
STORAGE**




47287

W : 36
H : 84
D : 23

- Upper/Lower Doors
- SA902 Molded Trays: 48

**TRAY
STORAGE**



47288

W : 48
H : 84
D : 23


- Upper/Lower Doors
- SA901 Molded Trays: 48



SA902: 10W x 3 1/2H x 19D

84H
TALL

**TRAY
STORAGE**



47281

W : 48
H : 84
D : 23

- SA902 Molded Trays: 64

**TRAY STORAGE/
SHELF**




47302

W : 36
H : 84
D : 23

- Upper/Lower Doors
- Adjustable Shelves: 3
- SA902 Molded Trays: 18


**TRAY
STORAGE**



47282


W : 48
H : 84
D : 23

- Doors
- SA902 Molded Trays: 64



Tray and bin racks have
captured runner rail to
reduce tray tipout.

**TRAY
STORAGE**



47289

W : 48
H : 84
D : 23

- Upper/Lower Doors
- SA902 Molded Trays: 64




SA910: 12 $\frac{1}{4}$ W x 3H x 16 $\frac{3}{4}$ D



SA911: 12 $\frac{1}{4}$ W x 6H x 16 $\frac{3}{4}$ D

84H
TALL

**BIN
STORAGE**




47242

W : 42
H : 84
D : 23

- SA910 Molded Translucent Bins: 48

**BIN
STORAGE**




47244

W : 42
H : 84
D : 23

- SA911 Molded Translucent Bins: 24

**BIN
STORAGE**




47243

W : 42
H : 84
D : 23

- Doors
- SA910 Molded Translucent Bins: 48


**BIN
STORAGE**



47245

W : 42
H : 84
D : 23

- Doors
- SA911 Molded Translucent Bins: 24



Tray and bin racks have
captured runner rail to
reduce tray tipout.

**FILE
DRAWER**

47326

W : 15, 18
H : 84
D : 23

- Adjustable Shelves: 2
- Full Extension File Drawers: 3
- Includes Crossbars for Side Filing
- Letter File Drawers (15W)
- Legal File Drawers (18W)

**FILE
DRAWER**

47322

W : 24, 30, 36
H : 84
D : 23


- Adjustable Shelves: 2
- Full Extension Lateral File Drawers: 3
- Includes Divider for Side Filing and File Bars for Front-to-Back Filing
- Drawer Box Depth: 440mm

**FILE
DRAWER**

47327

W : 15, 18
H : 84
D : 23

- Right Hinge Door
- Adjustable Shelves: 2
- Full Extension File Drawers: 3
- Includes Crossbars for Side Filing
- Letter File Drawers (15W)
- Legal File Drawers (18W)

**FILE
DRAWER**

47323

W : 24, 30, 36
H : 84
D : 23

- Doors
- Adjustable Shelves: 2
- Full Extension Lateral File Drawers: 3
- Includes Divider for Side Filing and File Bars for Front-to-Back Filing
- Drawer Box Depth: 440mm

**FILE
DRAWER**

47328

W : 15, 18
H : 84
D : 23

- Left Hinge Door
- Adjustable Shelves: 2
- Full Extension File Drawers: 3
- Includes Crossbars for Side Filing
- Letter File Drawers (15W)
- Legal File Drawers (18W)




**SMOCK
STORAGE**



47611

W : 18, 24
H : 84
D : 23

- Right Hinge Door
- Fixed Shelf
- SA364MN Double J Shelf Hooks: 14




**SMOCK
STORAGE**




47612

W : 18, 24
H : 84
D : 23

- Left Hinge Door
- Fixed Shelf
- SA364MN Double J Shelf Hooks: 14




**SMOCK
STORAGE**



47556

W : 36, 48
H : 84
D : 23

- Doors
- Center Divider
- Fixed Shelves
- SA364MN Double J Shelf Hooks: 28




**WARDROBE/
STORAGE**

47543

W : 15, 18, 21, 24
H : 84
D : 23

- Right Hinge Door
- Hanger Bar
- Adjustable Shelf (Factory Located)

**WARDROBE/
STORAGE**

47618

W : 36
H : 84
D : 23


- Locked Doors
- Full Height Divider
- Left: 5 Adjustable Shelves
- Right: Adjustable Shelf, Wardrobe Hook, Mirror

**WARDROBE/
STORAGE**

47544


W : 15, 18, 21, 24
H : 84
D : 23

- Left Hinge Door
- Hanger Bar
- Adjustable Shelf (Factory Located)

**WARDROBE/
STORAGE**

47603


W : 30, 36, 42, 48
H : 84
D : 23

- Locked Doors
- Full Height Divider
- Left: Adjustable Shelf, Wardrobe Hook, Mirror
- Right: 5 Adjustable Shelves

**WARDROBE/
STORAGE**

47617


W : 36
H : 84
D : 23

- Locked Doors
- Full Height Divider
- Left: Adjustable Shelf, Wardrobe Hook, Mirror
- Right: 5 Adjustable Shelves

**WARDROBE/
STORAGE**

47605


W : 30, 36, 42, 48
H : 84
D : 23

- Locked Doors
- Full Height Divider
- Left: 5 Adjustable Shelves
- Right: Adjustable Shelf, Wardrobe Hook, Mirror

**WARDROBE/
STORAGE**

47657


W : 30, 36
H : 84
D : 23

- Locked Doors
- Full Height Center Divider
- Left: Adjustable Shelf, Hanger Bar, Mirror
- Right: 3 Full Extension Locked File Drawers, Storage Compartment with 2 Adjustable Shelves
- Includes Crossbar for Side Filing
- Letter File Drawers (30W)
- Legal File Drawers (36W)

**WARDROBE/
STORAGE**

47607


W : 36, 42, 48
H : 84
D : 23

- Locked Doors
- Full Height Divider
- Left: 3 Adjustable Shelves, 2 Full Extension Letter File Drawers
- Right: Adjustable Shelf, Wardrobe Hook, Mirror

**WARDROBE/
STORAGE**

47653

W : 30, 36
H : 84
D : 23

- Locked Doors
- Full Height Center Divider
- Left: 3 Full Extension Locked File Drawers, Storage Compartment with 2 Adjustable Shelves
- Right: Adjustable Shelf, Hanger Bar, Mirror
- Includes Crossbar for Side Filing
- Letter File Drawers (30W)
- Legal File Drawers (36W)

**WARDROBE/
STORAGE**

47606

W : 36, 42, 48
H : 84
D : 23

- Locked Doors
- Full Height Divider
- Left: Adjustable Shelf, Wardrobe Hook, Mirror
- Right: 3 Adjustable Shelves, 2 Full Extension Letter File Drawers

**CHEMICAL
STORAGE**



47734

W : 15, 18
H : 84
D : 19, 23

- Right Hinge Door
- Black A-tech Shelves with Front Retaining Ledge: 5

**APPARATUS
ROD STORAGE**



47712

W : 15, 18
H : 84
D : 23

- Right Hinge Door
- Vertical Storage Compartments: 2
- Shelf
- Clamp Bin

**CHEMICAL
STORAGE**



47735

W : 15, 18
H : 84
D : 19, 23

- Left Hinge Door
- Black A-tech Shelves with Front Retaining Ledge: 5

**APPARATUS
ROD STORAGE**



47713

W : 15, 18
H : 84
D : 23

- Left Hinge Door
- Vertical Storage Compartments: 2
- Shelf
- Clamp Bin

**CHEMICAL
STORAGE**



47737

W : 30, 36
H : 84
D : 19, 23

- Doors
- Center Divider
- Black A-tech Shelves with Front Retaining Ledge: 10



**TORSO/SKELETON
STORAGE**

47722



W : 24
H : 84
D : 23

- Right Hinge Display Door
- Pullout Skeleton Rack
- Slideout Torso Shelf
- Adjustable Shelf


**TORSO/SKELETON
STORAGE**

47723



W : 24
H : 84
D : 23

- Left Hinge Display Door
- Pullout Skeleton Rack
- Slideout Torso Shelf
- Adjustable Shelf

**MICROSCOPE
STORAGE**

47728

W : 24
H : 84
D : 19, 23


- Right Hinge Door
- Black A-tech Shelving
- Vertical Dividers
- Compartments: 15

**MICROSCOPE
STORAGE**

47730

W : 24
H : 84
D : 19, 23

- Right Hinge Display Door
- Black A-tech Shelving
- Vertical Dividers
- Compartments: 15

**MICROSCOPE
STORAGE**

47729

W : 24
H : 84
D : 19, 23

- Left Hinge Door
- Black A-tech Shelving
- Vertical Dividers
- Compartments: 15

**MICROSCOPE
STORAGE**

47731

W : 24
H : 84
D : 19, 23


- Left Hinge Display Door
- Black A-tech Shelving
- Vertical Dividers
- Compartments: 15

**MICROSCOPE
STORAGE**

47732

W : 36
H : 84
D : 19, 23

- Doors
- Black A-tech Shelving
- Vertical Dividers
- Compartments: 20

**MICROSCOPE
STORAGE**

47733

W : 36
H : 84
D : 19, 23

- Display Doors
- Black A-tech Shelving
- Vertical Dividers
- Compartments: 20



SHELF



47101

W : 12, 15, 18, 21, 24, 30, 36
H : 90
D : 14, 17, 19, 23

- Fixed Shelf
- Adjustable Shelves: 4

SHELF



47129

W : 30, 36
H : 90
D : 14, 17, 19, 23

- Doors
- Fixed Shelf
- Adjustable Shelves: 4

SHELF




47120

W : 12, 15, 18, 21, 24
H : 90
D : 14, 17, 19, 23

- Right Hinge Door
- Fixed Shelf
- Adjustable Shelves: 4

SHELF




47126

W : 30, 36
H : 90
D : 14, 17, 19, 23

- Upper/Lower Doors
- Adjustable Shelves: 4

SHELF



47121

W : 12, 15, 18, 21, 24
H : 90
D : 14, 17, 19, 23

- Left Hinge Door
- Fixed Shelf
- Adjustable Shelves: 4

SHELF




47138

W : 42, 48
H : 90
D : 14, 17, 19, 23


- Center Divider
- Adjustable Shelves: 10

SHELF


47181

W : 42, 48
H : 90
D : 14, 17, 19, 23

- Doors
- Center Divider
- Adjustable Shelves: 10

SHELF/
DRAWER

47335

W : 30, 36, 42, 48
H : 90
D : 17, 19, 23


- Doors
- Center Divider
- Adjustable Shelves: 6
- Drawers: 3

SHELF


47149

W : 42, 48
H : 90
D : 17, 19, 23

- Sliding Doors
- Center Divider
- Adjustable Shelves: 10

SHELF/
DRAWER

47374

W : 30, 36, 42, 48
H : 90
D : 17, 19, 23

- Doors
- Center Divider
- Adjustable Shelves: 6
- Drawers: 4

SHELF


47150

W : 42, 48
H : 90
D : 17, 19, 23

- Upper/Lower Sliding Doors
- Center Dividers
- Adjustable Shelves: 8

DISPLAY




47130

W : 15, 18, 21, 24
H : 90
D : 17, 19, 23

- Right Hinge Display Door
- Fixed Shelf
- Adjustable Shelves: 4

DISPLAY



47132

W : 15, 18, 21, 24
H : 90
D : 17, 19, 23

- Right Hinge Display Door
- Fixed Shelf
- Adjustable Shelves: 4

DISPLAY




47131

W : 15, 18, 21, 24
H : 90
D : 17, 19, 23

- Left Hinge Display Door
- Fixed Shelf
- Adjustable Shelves: 4

DISPLAY



47133

W : 15, 18, 21, 24
H : 90
D : 17, 19, 23

- Left Hinge Display Door
- Fixed Shelf
- Adjustable Shelves: 4

DISPLAY



47137

W : 30, 36, 42, 48
H : 90
D : 17, 19, 23

- Display Doors
- Center Divider
- Adjustable Shelves: 10

DISPLAY



47140

W : 30, 36, 42, 48
H : 90
D : 17, 19, 23


- Display Doors
- Center Divider
- Adjustable Shelves: 10

DISPLAY

47182

W : 15, 18, 21, 24
H : 90
D : 17, 19, 23

- Right Hinge Upper/Lower Display Doors
- Adjustable Shelves: 4

DISPLAY

47135

W : 30, 36
H : 90
D : 17, 19, 23

- Upper Display Doors
- Lower Doors
- Adjustable Shelves: 4

DISPLAY

47183

W : 15, 18, 21, 24
H : 90
D : 17, 19, 23

- Left Hinge Upper/Lower Display Doors
- Adjustable Shelves: 4

DISPLAY

47154

W : 30, 36, 42, 48
H : 90
D : 17, 19, 23

- Upper Sliding Display Doors
- Lower Sliding Doors
- Center Dividers
- Adjustable Shelves: 8

DISPLAY

47184

W : 30, 36
H : 90
D : 17, 19, 23

- Upper/Lower Display Doors
- Adjustable Shelves: 4

DISPLAY

47158

W : 30, 36, 42, 48
H : 90
D : 17, 19, 23

- Upper Sliding Glass Doors
- Lower Doors
- Center Dividers
- Adjustable Shelves: 8

**DISPLAY/
DRAWER**

47339

W : 30, 36, 42, 48
H : 90
D : 17, 19, 23

- Display Doors
- Center Divider
- Adjustable Shelves: 6
- Drawers: 3

**DISPLAY/
DRAWER**

47379

W : 30, 36, 42, 48
H : 90
D : 17, 19, 23

- Display Doors
- Center Divider
- Adjustable Shelves: 6
- Drawers: 4

**DISPLAY/
DRAWER**

47342

W : 30, 36, 42, 48
H : 90
D : 17, 19, 23

- Sliding Display Doors
- Center Divider
- Adjustable Shelves: 6
- Drawers: 3

**DISPLAY/
DRAWER**

47378

W : 30, 36, 42, 48
H : 90
D : 17, 19, 23

- Sliding Display Doors
- Center Divider
- Adjustable Shelves: 6
- Drawers: 4

**DISPLAY/
DRAWER**

47343

W : 30, 36, 42, 48
H : 90
D : 17, 19, 23

- Sliding Glass Doors
- Center Divider
- Adjustable Shelves: 6
- Drawers: 3

**DISPLAY/
DRAWER**

47383

W : 30, 36, 42, 48
H : 90
D : 17, 19, 23

- Sliding Glass Doors
- Center Divider
- Adjustable Shelves: 6
- Drawers: 4

**EXPANDED
METAL**

47106

W : 15, 18, 21, 24
H : 90
D : 17, 19, 23

- Right Hinge Framed Expanded Metal Door
- Adjustable Shelves: 5
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**

47109

W : 30, 36, 42, 48
H : 90
D : 17, 19, 23

- Framed Expanded Metal Doors
- Center Divider
- Adjustable Shelves: 10
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**

47107

W : 15, 18, 21, 24
H : 90
D : 17, 19, 23

- Left Hinge Framed Expanded Metal Door
- Adjustable Shelves: 5
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**


47108

W : 30, 36
H : 90
D : 17, 19, 23

- Framed Expanded Metal Doors
- Adjustable Shelves: 5
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)



CUBICLE



47213

W : 21, 24
H : 90
D : 17, 19, 23

- 3/4" Adjustable Cubicle Shelves: 14
- Compartments: 16
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE



47216

W : 30, 36, 42, 48
H : 90
D : 17, 19, 23

- 3/4" Adjustable Cubicle Shelves: 21
- Compartments: 24
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W), 13" (42W), 15" (48W)

CUBICLE



47214

W : 21, 24
H : 90
D : 17, 19, 23

- Right Hinge Door
- 3/4" Adjustable Cubicle Shelves: 14
- Compartments: 16
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE




47217

W : 30, 36, 42, 48
H : 90
D : 17, 19, 23

- Doors
- 3/4" Adjustable Cubicle Shelves: 21
- Compartments: 24
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W), 13" (42W), 15" (48W)

CUBICLE



47215

W : 21, 24
H : 90
D : 17, 19, 23

- Left Hinge Door
- 3/4" Adjustable Cubicle Shelves: 14
- Compartments: 16
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)


CUBICLE



47218

W : 36, 42, 48
H : 90
D : 17, 19, 23

- 3/4" Adjustable Cubicle Shelves: 28
- Compartments: 32
- Compartment Spacing: 3"H min.
- Compartment Widths: 8" (36W), 9 1/2" (42W), 11" (48W)

CUBICLE

47219

W : 36, 42, 48
H : 90
D : 17, 19, 23

- Doors
- 3/4" Adjustable Cubicle Shelves: 28
- Compartments: 32
- Compartment Spacing: 3"H min.
- Compartment Widths: 8" (36W), 9 1/2" (42W), 11" (48W)

**CUBICLE/
SHELF**

47247

W : 36
H : 90
D : 17, 19, 23

- Doors
- Adjustable Shelves: 3
- 3/4" Adjustable Cubicle Shelves: 12
- Compartments: 16
- Compartment Spacing: 3"H min.
- Compartment Width: 8"

**CUBICLE/
SHELF**

47232

W : 30, 36
H : 90
D : 17, 19, 23


- Doors
- Adjustable Shelves: 3
- 3/4" Adjustable Cubicle Shelves: 9
- Compartments: 12
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W)

**CUBICLE/
SHELF**

47246

W : 36
H : 90
D : 17, 19, 23

- Upper/Lower Doors
- Adjustable Shelves: 3
- 3/4" Adjustable Cubicle Shelves: 12
- Compartments: 16
- Compartment Spacing: 3"H min.
- Compartment Width: 8"


**CUBICLE/
SHELF**

47231

W : 30, 36
H : 90
D : 17, 19, 23

- Upper/Lower Doors
- Adjustable Shelves: 3
- 3/4" Adjustable Cubicle Shelves: 9
- Compartments: 12
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W)



**BLIND
CORNER**



47200

W : 33, 36, 39
H : 90
D : 17

W : 36, 39, 42
H : 90
D : 19

W : 39, 42, 45, 48
H : 90
D : 23

- Right Hinge Door
- Fixed Shelf
- Adjustable Shelves: 4


**SMOCK
STORAGE**




47611

W : 18, 24
H : 90
D : 23

- Right Hinge Door
- Fixed Shelf
- SA364MN Double J Shelf Hooks: 14



**BLIND
CORNER**



47201

W : 33, 36, 39
H : 90
D : 17

W : 36, 39, 42
H : 90
D : 19

W : 39, 42, 45, 48
H : 90
D : 23

- Left Hinge Door
- Fixed Shelf
- Adjustable Shelves: 4


**SMOCK
STORAGE**



47612

W : 18, 24
H : 90
D : 23

- Left Hinge Door
- Fixed Shelf
- SA364MN Double J Shelf Hooks: 14




**SMOCK
STORAGE**



47556

W : 36, 48
H : 90
D : 23

- Doors
- Center Divider
- Fixed Shelves: 2
- SA364MN Double J Shelf Hooks: 28




**WARDROBE/
STORAGE**

47543

W : 15, 18, 21, 24
H : 90
D : 23

- Right Hinge Door
- Hanger Bar
- Adjustable Shelf (Factory Located)

**WARDROBE/
STORAGE**

47618

W : 36
H : 90
D : 23


- Locked Doors
- Full Height Divider
- Left: 5 Adjustable Shelves
- Right: Adjustable Shelf, Wardrobe Hook, Mirror

**WARDROBE/
STORAGE**

47544


W : 15, 18, 21, 24
H : 90
D : 23

- Left Hinge Door
- Hanger Bar
- Adjustable Shelf (Factory Located)

**WARDROBE/
STORAGE**

47603


W : 30, 36, 42, 48
H : 90
D : 23

- Locked Doors
- Full Height Divider
- Left: Adjustable Shelf, Wardrobe Hook, Mirror
- Right: 5 Adjustable Shelves

**WARDROBE/
STORAGE**

47617


W : 36
H : 90
D : 23

- Locked Doors
- Full Height Divider
- Left: Adjustable Shelf, Wardrobe Hook, Mirror
- Right: 5 Adjustable Shelves

**WARDROBE/
STORAGE**

47605


W : 30, 36, 42, 48
H : 90
D : 23

- Locked Doors
- Full Height Divider
- Left: 5 Adjustable Shelves
- Right: Adjustable Shelf, Wardrobe Hook, Mirror

**WARDROBE/
STORAGE**

47657


W : 30, 36
H : 90
D : 23

- Locked Doors
- Full Height Center Divider
- Left: Adjustable Shelf, Hanger Bar, Mirror
- Right: 3 Full Extension Locked File Drawers, Storage Compartment with 2 Adjustable Shelves
- Includes Crossbar for Side Filing
- Letter File Drawers (30W)
- Legal File Drawers (36W)

**WARDROBE/
STORAGE**

47607


W : 36, 42, 48
H : 90
D : 23

- Locked Doors
- Full Height Divider
- Left: 3 Adjustable Shelves, 2 Full Extension Letter File Drawers
- Right: Adjustable Shelf, Wardrobe Hook, Mirror

**WARDROBE/
STORAGE**

47653

W : 30, 36
H : 90
D : 23

- Locked Doors
- Full Height Center Divider
- Left: 3 Full Extension Locked File Drawers, Storage Compartment with 2 Adjustable Shelves
- Right: Adjustable Shelf, Hanger Bar, Mirror
- Includes Crossbar for Side Filing
- Letter File Drawers (30W)
- Legal File Drawers (36W)

**WARDROBE/
STORAGE**

47606

W : 36, 42, 48
H : 90
D : 23


- Locked Doors
- Full Height Divider
- Left: Adjustable Shelf, Wardrobe Hook, Mirror
- Right: 3 Adjustable Shelves, 2 Full Extension Letter File Drawers

**APPARATUS
ROD STORAGE**

47712

W : 15, 18
H : 90
D : 23

- Right Hinge Door
- Vertical Storage Compartments: 2
- Shelf
- Clamp Bin

**ADA BARRIER FREE
SAFETY STATION**

47699

W : 36
H : 90
D : 24


- 2 Fixed Shelves
- Solid Surface Top and Apron
- Front Access Panel
- Shower Head
- 4-Head Eye/Face Wash
- Push-Down Valve
- Recessed Stainless Steel Bowl

**APPARATUS
ROD STORAGE**

47713

W : 15, 18
H : 90
D : 23

- Left Hinge Door
- Vertical Storage Compartments: 2
- Shelf
- Clamp Bin

**SAFETY
STATION**

47700

W : 36
H : 90
D : 23

- Fixed Shelf
- Solid Surface Top
- Front Access Panel
- Shower Head
- Back-Mounted 2-Head Eye/Face Wash
- ABS Bowl

**MICROSCOPE
STORAGE**




47728

W : 24
H : 90
D : 19, 23

- Right Hinge Door
- Black A-tech Shelving
- Vertical Dividers
- Compartments: 15

**MICROSCOPE
STORAGE**




47730

W : 24
H : 90
D : 19, 23

- Right Hinge Display Door
- Black A-tech Shelving
- Vertical Dividers
- Compartments: 15

**MICROSCOPE
STORAGE**




47729

W : 24
H : 90
D : 19, 23

- Left Hinge Door
- Black A-tech Shelving
- Vertical Dividers
- Compartments: 15

**MICROSCOPE
STORAGE**



47731

W : 24
H : 90
D : 19, 23

- Left Hinge Display Door
- Black A-tech Shelving
- Vertical Dividers
- Compartments: 15

**MICROSCOPE
STORAGE**



47732

W : 36
H : 90
D : 19, 23

- Doors
- Black A-tech Shelving
- Vertical Dividers
- Compartments: 20

**MICROSCOPE
STORAGE**



47733


W : 36
H : 90
D : 19, 23

- Display Doors
- Black A-tech Shelving
- Vertical Dividers
- Compartments: 20





Stack-On units are designed to be placed on top of standard tall cabinets, allowing additional upper storage utilization.

Advantage Scientific cabinets can be produced to a maximum 90H nominal in a one-piece vertical assembly. When taller cabinet assemblies are desired, Stack-On units provide the solution. Stack-Ons may be added to 84H nominal (2137mm / 84.13") and 90H nominal (2297mm / 90.43") Tall cabinets.




Stack-On units may also be utilized to span more than one Tall cabinet below. For example, two 24"-wide Tall units could use a single 48"-wide Stack-On. Note that, due to side hinge reveals, vertical gaps will align, but width spacing of gaps may vary due to hinge hardware selected.




STACK-ON WITH 84H TALL



STACK-ON WITH 90H TALL



SHELF




43500

W : 15, 18, 21, 24, 30, 36
H : 17
D : 14, 17, 19, 23

- Open Storage
- For Use with 84H and 90H Tall Cabinets

SHELF



43508

W : 42, 48
H : 17
D : 14, 17, 19, 23

- Center Divider
- For Use with 84H and 90H Tall Cabinets

SHELF



43502

W : 15, 18, 21, 24
H : 17
D : 14, 17, 19, 23

- Right Hinge Door
- For Use with 84H and 90H Tall Cabinets

SHELF



43510

W : 42, 48
H : 17
D : 14, 17, 19, 23

- Doors
- Center Divider
- For Use with 84H and 90H Tall Cabinets

SHELF



43504

W : 15, 18, 21, 24
H : 17
D : 14, 17, 19, 23

- Left Hinge Door
- For Use with 84H and 90H Tall Cabinets

SHELF



43520

W : 42, 48
H : 17
D : 17, 19, 23

- Sliding Doors
- Center Divider
- For Use with 84H and 90H Tall Cabinets

SHELF



43506

W : 30, 36
H : 17
D : 14, 17, 19, 23

- Doors
- For Use with 84H and 90H Tall Cabinets

DISPLAY



43524

W : 42, 48
H : 17
D : 17, 19, 23

- Sliding Glass Doors
- Center Divider
- For Use with 84H and 90H Tall Cabinets

SHELF




43500

W : 15, 18, 21, 24, 30, 36
H : 23
D : 14, 17, 19, 23

- Adjustable Shelf
- For Use with 84H and 90H Tall Cabinets

SHELF



43508

W : 42, 48
H : 23
D : 14, 17, 19, 23

- Center Divider
- Adjustable Shelves: 2
- For Use with 84H and 90H Tall Cabinets

SHELF




43502

W : 15, 18, 21, 24
H : 23
D : 14, 17, 19, 23

- Right Hinge Door
- Adjustable Shelf
- For Use with 84H and 90H Tall Cabinets

SHELF



43510

W : 42, 48
H : 23
D : 14, 17, 19, 23

- Doors
- Center Divider
- Adjustable Shelves: 2
- For Use with 84H and 90H Tall Cabinets

SHELF



43504

W : 15, 18, 21, 24
H : 23
D : 14, 17, 19, 23

- Left Hinge Door
- Adjustable Shelf
- For Use with 84H and 90H Tall Cabinets

SHELF



43520

W : 42, 48
H : 23
D : 17, 19, 23

- Sliding Doors
- Center Divider
- Adjustable Shelves: 2
- For Use with 84H and 90H Tall Cabinets

SHELF



43506

W : 30, 36
H : 23
D : 14, 17, 19, 23

- Doors
- Adjustable Shelf
- For Use with 84H and 90H Tall Cabinets



DISPLAY



43512

W : 15, 18, 21, 24
H : 23
D : 14, 17, 19, 23

- Right Hinge Display Door
- Adjustable Shelf
- For Use with 84H and 90H Tall Cabinets

DISPLAY



43522

W : 42, 48
H : 23
D : 17, 19, 23

- Sliding Display Doors
- Center Divider
- Adjustable Shelves: 2
- For Use with 84H and 90H Tall Cabinets

DISPLAY




43514

W : 15, 18, 21, 24
H : 23
D : 14, 17, 19, 23

- Left Hinge Display Door
- Adjustable Shelf
- For Use with 84H and 90H Tall Cabinets

DISPLAY



43524

W : 42, 48
H : 23
D : 17, 19, 23

- Sliding Glass Doors
- Center Divider
- Adjustable Shelves: 2
- For Use with 84H and 90H Tall Cabinets

DISPLAY



43518

W : 30, 36
H : 23
D : 14, 17, 19, 23

- Display Doors
- Adjustable Shelf
- For Use with 84H and 90H Tall Cabinets

DISPLAY



43516

W : 42, 48
H : 23
D : 14, 17, 19, 23

- Display Doors
- Center Divider
- Adjustable Shelves: 2
- For Use with 84H and 90H Tall Cabinets

SHELF




43500

W : 15, 18, 21, 24, 30, 36
H : 29
D : 14, 17, 19, 23

- Adjustable Shelves: 2
- For Use with 84H and 90H Tall Cabinets

SHELF



43508

W : 42, 48
H : 29
D : 14, 17, 19, 23

- Center Divider
- Adjustable Shelves: 4
- For Use with 84H and 90H Tall Cabinets

SHELF




43502

W : 15, 18, 21, 24
H : 29
D : 14, 17, 19, 23

- Right Hinge Door
- Adjustable Shelves: 2
- For Use with 84H and 90H Tall Cabinets

SHELF



43510

W : 42, 48
H : 29
D : 14, 17, 19, 23

- Doors
- Center Divider
- Adjustable Shelves: 4
- For Use with 84H and 90H Tall Cabinets

SHELF




43504

W : 15, 18, 21, 24
H : 29
D : 14, 17, 19, 23

- Left Hinge Door
- Adjustable Shelves: 2
- For Use with 84H and 90H Tall Cabinets

SHELF




43520

W : 42, 48
H : 29
D : 17, 19, 23

- Sliding Doors
- Center Divider
- Adjustable Shelves: 4
- For Use with 84H and 90H Tall Cabinets

SHELF



43506

W : 30, 36
H : 29
D : 14, 17, 19, 23

- Doors
- Adjustable Shelves: 2
- For Use with 84H and 90H Tall Cabinets



DISPLAY




43512

W : 15, 18, 21, 24
H : 29
D : 14, 17, 19, 23

- Right Hinge Display Door
- Adjustable Shelves: 2
- For Use with 84H and 90H Tall Cabinets

DISPLAY



43522

W : 42, 48
H : 29
D : 17, 19, 23

- Sliding Display Doors
- Center Divider
- Adjustable Shelves: 4
- For Use with 84H and 90H Tall Cabinets

DISPLAY




43514

W : 15, 18, 21, 24
H : 29
D : 14, 17, 19, 23

- Left Hinge Display Door
- Adjustable Shelves: 2
- For Use with 84H and 90H Tall Cabinets

DISPLAY




43524

W : 42, 48
H : 29
D : 17, 19, 23

- Sliding Glass Doors
- Center Divider
- Adjustable Shelves: 4
- For Use with 84H and 90H Tall Cabinets

DISPLAY



43518

W : 30, 36
H : 29
D : 14, 17, 19, 23

- Display Doors
- Adjustable Shelves: 2
- For Use with 84H and 90H Tall Cabinets

DISPLAY




43516

W : 42, 48
H : 29
D : 14, 17, 19, 23

- Display Doors
- Center Divider
- Adjustable Shelves: 4
- For Use with 84H and 90H Tall Cabinets

SHELF




43500

W : 15, 18, 21, 24, 30, 36
H : 35
D : 14, 17, 19, 23

- Adjustable Shelves: 2
- For Use with 84H Tall Cabinets

SHELF



43508

W : 42, 48
H : 35
D : 14, 17, 19, 23

- Center Divider
- Adjustable Shelves: 2
- For Use with 84H Tall Cabinets

SHELF



43502

W : 15, 18, 21, 24
H : 35
D : 14, 17, 19, 23

- Right Hinge Door
- Adjustable Shelves: 2
- For Use with 84H Tall Cabinets

SHELF



43510

W : 42, 48
H : 35
D : 14, 17, 19, 23

- Doors
- Center Divider
- Adjustable Shelves: 4
- For Use with 84H Tall Cabinets

SHELF




43504

W : 15, 18, 21, 24
H : 35
D : 14, 17, 19, 23

- Left Hinge Door
- Adjustable Shelves: 2
- For Use with 84H Tall Cabinets

SHELF




43520

W : 42, 48
H : 35
D : 17, 19, 23

- Sliding Doors
- Center Divider
- Adjustable Shelves: 4
- For Use with 84H Tall Cabinets

SHELF



43506

W : 30, 36
H : 35
D : 14, 17, 19, 23

- Doors
- Adjustable Shelves: 2
- For Use with 84H Tall Cabinets



DISPLAY



43512

W : 15, 18, 21, 24
H : 35
D : 14, 17, 19, 23

- Right Hinge Display Door
- Adjustable Shelves: 2
- For Use with 84H Tall Cabinets

DISPLAY



43522

W : 42, 48
H : 35
D : 17, 19, 23

- Sliding Display Doors
- Center Divider
- Adjustable Shelves: 4
- For Use with 84H Tall Cabinets

DISPLAY




43514

W : 15, 18, 21, 24
H : 35
D : 14, 17, 19, 23

- Left Hinge Display Door
- Adjustable Shelves: 2
- For Use with 84H Tall Cabinets

DISPLAY



43524

W : 42, 48
H : 35
D : 17, 19, 23

- Sliding Glass Doors
- Center Divider
- Adjustable Shelves: 4
- For Use with 84H Tall Cabinets

DISPLAY



43518

W : 30, 36
H : 35
D : 14, 17, 19, 23

- Display Doors
- Adjustable Shelves: 2
- For Use with 84H Tall Cabinets

DISPLAY



43516

W : 42, 48
H : 35
D : 14, 17, 19, 23

- Display Doors
- Center Divider
- Adjustable Shelves: 4
- For Use with 84H Tall Cabinets

SOFFIT TRIM



25801

W : 96
H : 00
D : 03

- Field Trim to Fit
- Actual Thickness: 3/4"
- Actual Depth: 3"
- Edgebanding: One Long Edge in Casebody Trim Color and 0.5mm Thickness

VERTICAL SCRIBE




25802

W : 00
H : 96
D : 14, 17, 19, 23

- Field Trim to Fit
- Actual Thickness: 3/4"
- Actual Depth: 3"
- Edgebanding: One Long Edge in Casebody Trim Color and Thickness

FRONT FILLER




25803

W : 02, 03, 06, 09, 12
H : 84, 90
D : 14, 17, 19, 23

- Field Trim to Fit
- Edgebanding: None

FRONT / TOP FILLER





25811

W : 02, 03, 06, 09, 12
H : 84, 90
D : 14, 17, 19, 23

- Field Trim to Fit
- Edgebanding: None

DIAGONAL CORNER FILLERS



25804 RIGHT


25806 LEFT

W : 00
H : 84, 90
D : 14, 17, 19, 23

- Field Trim to Fit
- Depths Listed are Depths of Cabinets to which Filler is Attached
- Edgebanding: One Long Edge in Casebody Trim Color and 0.5mm Thickness



CORNER FILLER




25805

W : 03, 06, 09, 12
H : 84, 90
D : 14, 17, 19, 23

- Field Trim to Fit
- Depths Listed are Depths of Cabinets to which Filler is Attached
- Edgebanding: None

SIDE FILLER




25807

W : 00
H : 84, 90
D : 04, 06, 08, 12, 15

- Field Trim to Fit
- Edgebanding: One Long Edge in Casebody Trim Color and Thickness

TOP CORNER CLOSURE




25809

W : 16, 18, 21, 24, 26, 29, 31, 40, 42
H : 00
D : 16, 18, 21, 24, 26, 29, 31, 40, 42

- Field Trim to Fit
- Edgebanding: One Short and One Long Edge in Casebody Trim Color and Thickness

TOP BACK CLOSURE




25810

W : 12, 15, 18, 21, 24, 30, 36, 42, 48
H : 00
D : 04, 06, 08, 12, 15

- Horizontal Grain
- Field Trim to Fit
- Edgebanding: None

Shown with
25807 Side
Filler

END SOFFIT CLOSURE




25820

W : 00
H : 12, 18, 24
D : 14, 17, 19, 23

- Vertical Grain
- Edgebanding: All Edges in Casebody Trim Color and 0.5mm Thickness

SOFFIT EXTENSION




25824

W : 12, 15, 18, 21, 24, 30, 36, 42, 48
H : 6"-34" (1" increments)
D : 14, 17, 19, 23

- Edgebanding (Box): All Edges in Casebody Trim Color and Thickness
- Edgebanding (Soffit Panel): All Edges in Door/Drawer Trim Color and Thickness

SOFFIT FASCIA




25822

W : 12, 15, 18, 21, 24, 30, 36, 42, 48, 96
H : 12, 18, 24
D : 00

- Horizontal Grain
- Edgebanding: One Long and Two Short Edges in Casebody Trim Color and 0.5mm Thickness

INTERMEDIATE SOFFIT BRACKET




25821

W : 04
H : 12
D : 12

- Provides Added Support for Longer 25822 and 25823 Soffit Fascias
- Edgebanding: None

SOFFIT FASCIA




25823

W : 12, 15, 18, 21, 24, 30, 36, 42, 48
H : 12, 18, 24
D : 00

- Vertical Grain
- Edgebanding: One Short and Two Long Edges in Casebody Trim Color and 0.5mm Thickness



OVERLAY TOP



25830


W : 12, 15, 18, 21, 24, 30,
36, 42, 48, 96

H : 00


D : 14, 17, 19, 23

- Horizontal Grain
- Shipped in Lengths for Field Trimming
- Edgebanding: One Long and Two Short Edges in Casebody Trim Color and Thickness

DIAGONAL CORNER TOP



Shown with 25804/25806
Diagonal Corner Fillers



25812


W : 24, 30, 36, 42, 48

H : 00

D : 14, 17, 19, 23

- Horizontal Grain
- Edgebanding: One Short and One Long Edge in Casebody Trim Color and Thickness

OVERLAY BACK FULL HEIGHT



25813


W : 12, 15, 18, 21, 24, 30,
36, 42, 48

H : 84, 90

D : 00

- Vertical Grain
- Full Height
- Order by Cabinet Height
- Edgebanding: All Edges in Casebody Trim Color and Thickness

OVERLAY BACK ABOVE KICK



25814


W : 12, 15, 18, 21, 24, 30,
36, 42, 48

H : 84, 90

D : 00

- Vertical Grain
- Above Kick
- Bottom Edge Aligns with Toe Kick
- Order by Cabinet Height
- Edgebanding: All Edges in Casebody Trim Color and Thickness

FILE CROSSBARS




10896

W : 15
H : 01
D : 00

- Pair of Crossbars for 15W File Drawers (letter)
- Allow Side Filing

FILE CROSSBARS




10897

W : 18
H : 01
D : 00

- Pair of Crossbars for 18W File Drawers (legal)
- Allow Side Filing

FILE BARS



10940

W : 00
H : 01
D : 18

- 1 Pair of File Bars for Lateral File Drawers
- Allow Front-to-Back Filing
- 1 Pair Required for 24W-27W Lateral File Drawers
- 2 Pairs Required for 30W-42W Lateral File Drawers