

18H	3	SHELF DISPLAY	3 3
25H	4	SHELF DISPLAY EXPANDED METAL CUBICLE CORNER DISPENSING	4 5 6 7 8 9
30H	10	SHELF DISPLAY EXPANDED METAL CUBICLE CORNER DISPENSING	10 11 12 13 14 15
36H	16	SHELF MICROWAVE SHELF DISPLAY EXPANDED METAL CUBICLE CORNER DISPENSING	16 16 17 18 19 20 21
42H	22	SHELF DISPLAY EXPANDED METAL CUBICLE CORNER	22 23 24 25 26
47H	27	SHELF DISPLAY EXPANDED METAL CUBICLE CORNER	27 28 29 30 31
INSTALLATION ACCESSORIES	32	SOFFIT TRIM VERTICAL SCRIBE FILLERS CLOSURES SOFFIT FASCIAS SOFFIT EXTENSION VALANCE RAIL	32 32 32 32 32 32 33
DESIGN DETAILS	34		



**ADVANTAGE
SCIENTIFIC**
BY STEVENS

VALANCE

SHELF


45601

W : 12, 14, 15, 17, 18, 20, 21, 23, 24, 28,
30, 34, 36

H : 18

D : 14, 17

- Open Storage

SHELF

45638

W : 24, 28, 30, 34, 36, 40, 42, 46, 48

H : 18

D : 14, 17

- Center Divider

SHELF

45620

W : 12, 14, 15, 17, 18

H : 18

D : 14, 17

- Right Hinge Door

SHELF

45681

W : 24, 28, 30, 34, 36

H : 18

D : 14, 17

- Doors
- Center Divider

SHELF

45621

W : 12, 14, 15, 17, 18

H : 18

D : 14, 17

- Left Hinge Door

SHELF

45650

W : 36, 40, 42, 46, 48

H : 18

D : 14, 17

- Sliding Doors
- Center Divider

SHELF

45629

W : 24, 28, 30, 34, 36

H : 18

D : 14, 17

- Doors

DISPLAY

45658

W : 36, 40, 42, 46, 48


H : 18

D : 14, 17

- Sliding Glass Doors
- Center Divider



SHELF




45601

W : 12, 14, 15, 17, 18, 20, 21, 23, 24, 28,
30, 34, 36
H : 25
D : 14, 17

- Open Storage
- Adjustable Shelf

SHELF



45638

W : 24, 28, 30, 34, 36, 40, 42, 46, 48
H : 25
D : 14, 17

- Center Divider
- Adjustable Shelves: 2

SHELF



45620

W : 12, 14, 15, 17, 18, 20, 21, 23, 24
H : 25
D : 14, 17

- Right Hinge Door
- Adjustable Shelf

SHELF



45681

W : 40, 42, 46, 48
H : 25
D : 14, 17

- Doors
- Center Divider
- Adjustable Shelves: 2

SHELF



45621

W : 12, 14, 15, 17, 18, 20, 21, 23, 24
H : 25
D : 14, 17

- Left Hinge Door
- Adjustable Shelf

SHELF



45650

W : 36, 40, 42, 46, 48
H : 25
D : 14, 17

- Sliding Doors
- Center Divider
- Adjustable Shelves: 2

SHELF



45629

W : 24, 28, 30, 34, 36
H : 25
D : 14, 17


- Doors
- Adjustable Shelf

DISPLAY

45630

W : 15, 17, 18, 20, 21, 23, 24
H : 25
D : 14, 17

- Right Hinge Display Door
- Adjustable Shelf

DISPLAY

45654

W : 36, 40, 42, 46, 48
H : 25
D : 14, 17


- Sliding Display Doors
- Center Divider
- Adjustable Shelves: 2

DISPLAY

45631

W : 15, 17, 18, 20, 21, 23, 24
H : 25
D : 14, 17

- Left Hinge Display Door
- Adjustable Shelf

DISPLAY

45658

W : 36, 40, 42, 46, 48
H : 25
D : 14, 17

- Sliding Glass Doors
- Center Divider
- Adjustable Shelves: 2

DISPLAY

45636

W : 30, 34, 36
H : 25
D : 14, 17

- Display Doors
- Adjustable Shelf

DISPLAY

45682

W : 40, 42, 46, 48
H : 25
D : 14, 17

- Display Doors
- Center Divider
- Adjustable Shelves: 2



**EXPANDED
METAL**



45606

W : 15, 17, 18, 20, 21, 23, 24
H : 25
D : 14, 17

- Right Hinge Framed Expanded Metal Door
- Adjustable Shelf
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**



45607

W : 15, 17, 18, 20, 21, 23, 24
H : 25
D : 14, 17

- Left Hinge Framed Expanded Metal Door
- Adjustable Shelf
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**



45608

W : 30, 34, 36
H : 25
D : 14, 17

- Framed Expanded Metal Doors
- Adjustable Shelf
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**



45609

W : 40, 42, 46, 48
H : 25
D : 14, 17

- Framed Expanded Metal Doors
- Center Divider
- Adjustable Shelves: 2
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

CUBICLE

45720

W : 21, 24
H : 25
D : 14, 17

- 3/4" Adjustable Cubicle Shelves: 4
- Compartments: 6
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE

45722

W : 21, 24
H : 25
D : 14, 17


- Right Hinge Door
- 3/4" Adjustable Cubicle Shelves: 4
- Compartments: 6
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE

45723

W : 21, 24
H : 25
D : 14, 17

- Left Hinge Door
- 3/4" Adjustable Cubicle Shelves: 4
- Compartments: 6
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE

45726

W : 30, 36, 40, 42, 46, 48
H : 25
D : 14, 17


- 3/4" Adjustable Cubicle Shelves: 6
- Compartments: 9
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W), 12 1/4" (40W), 13" (42W), 14 1/4" (46W), 15" (48W)

CUBICLE

45727

W : 30, 36, 40, 42, 46, 48
H : 25
D : 14, 17

- Doors
- 3/4" Adjustable Cubicle Shelves: 6
- Compartments: 9
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W), 12 1/4" (40W), 13" (42W), 14 1/4" (46W), 15" (48W)

CUBICLE

45739

W : 36, 40, 42, 46, 48
H : 25
D : 14, 17

- 3/4" Adjustable Cubicle Shelves: 8
- Compartments: 12
- Compartment Spacing: 3"H min.
- Compartment Widths: 8" (36W), 9" (40W), 9 1/2" (42W), 10 1/2" (46W), 11" (48W)

CUBICLE


45740

W : 36, 40, 42, 46, 48
H : 25
D : 14, 17

- Doors
- 3/4" Adjustable Cubicle Shelves: 8
- Compartments: 12
- Compartment Spacing: 3"H min.
- Compartment Widths: 8" (36W), 9" (40W), 9 1/2" (42W), 10 1/2" (46W), 11" (48W)



**BLIND
CORNER**




45700

W : 33, 36, 39
H : 25
D : 14

W : 36, 39, 42
H : 25
D : 17

- Right Hinge Door
- Adjustable Shelf

**DIAGONAL
CORNER**




45708

W : 25
H : 25
D : 14

W : 28
H : 25
D : 17

- Fixed Shelf

**BLIND
CORNER**



45701

W : 33, 36, 39
H : 25
D : 14

W : 36, 39, 42
H : 25
D : 17

- Left Hinge Door
- Adjustable Shelf

**DIAGONAL
CORNER**




45711

W : 25
H : 25
D : 14

W : 28
H : 25
D : 17

- Right Hinge Door
- Fixed Shelf

**BLIND
CORNER**



45702

W : 33, 36, 39
H : 25
D : 14

W : 36, 39, 42
H : 25
D : 17

- Right Hinge Display Door
- Adjustable Shelf

**DIAGONAL
CORNER**




45712

W : 25
H : 25
D : 14

W : 28
H : 25
D : 17

- Left Hinge Door
- Fixed Shelf

**BLIND
CORNER**



45703

W : 33, 36, 39
H : 25
D : 14

W : 36, 39, 42
H : 25
D : 17

- Left Hinge Display Door
- Adjustable Shelf

DISPENSING


45622

W : 21, 24
H : 25
D : 14, 17

- Right Hinge Door
- Towel Dispenser
- Adjustable Shelf

DISPENSING


45624

W : 18, 21, 24
H : 25
D : 14, 17

- Right Hinge Door
- Towel and Glove/Tissue Dispensers


DISPENSING


45623

W : 21, 24
H : 25
D : 14, 17

- Left Hinge Door
- Towel Dispenser
- Adjustable Shelf


DISPENSING


45625

W : 18, 21, 24
H : 25
D : 14, 17

- Left Hinge Door
- Towel and Glove/Tissue Dispensers


DISPENSING


45626

W : 36
H : 25
D : 14, 17


- Doors
- Towel Dispenser
- Adjustable Shelf

DISPENSING



45627

W : 36
H : 25
D : 14, 17

- Doors
- Towel and Glove/Tissue Dispensers




TOWEL DISPENSER
Bottom-access dispensing of paper towels.



GLOVE/TOISSUE DISPENSER
Face-of-door access for gloves or tissues.



SHELF




45601

W : 12, 14, 15, 17, 18, 20, 21, 23, 24, 28,
30, 34, 36
H : 30
D : 14, 17

- Open Storage
- Adjustable Shelf

SHELF



45638

W : 24, 28, 30, 34, 36, 40, 42, 46, 48
H : 30
D : 14, 17

- Center Divider
- Adjustable Shelves: 2

SHELF




45620

W : 12, 14, 15, 17, 18, 20, 21, 23, 24
H : 30
D : 14, 17

- Right Hinge Door
- Adjustable Shelf

SHELF



45681

W : 40, 42, 46, 48
H : 30
D : 14, 17

- Doors
- Center Divider
- Adjustable Shelves: 2

SHELF




45621

W : 12, 14, 15, 17, 18, 20, 21, 23, 24
H : 30
D : 14, 17

- Left Hinge Door
- Adjustable Shelf

SHELF



45650

W : 36, 40, 42, 46, 48
H : 30
D : 14, 17

- Sliding Doors
- Center Divider
- Adjustable Shelves: 2

SHELF



45629

W : 24, 28, 30, 34, 36
H : 30
D : 14, 17

- Doors
- Adjustable Shelf

DISPLAY

45630

W : 15, 17, 18, 20, 21, 23, 24
H : 30
D : 14, 17

- Right Hinge Display Door
- Adjustable Shelf

DISPLAY

45654

W : 36, 40, 42, 46, 48
H : 30
D : 14, 17

- Sliding Display Doors
- Center Divider
- Adjustable Shelves: 2

DISPLAY

45631

W : 15, 17, 18, 20, 21, 23, 24
H : 30
D : 14, 17

- Left Hinge Display Door
- Adjustable Shelf

DISPLAY

45658

W : 36, 40, 42, 46, 48
H : 30
D : 14, 17


- Sliding Glass Doors
- Center Divider
- Adjustable Shelves: 2

DISPLAY

45636

W : 30, 34, 36
H : 30
D : 14, 17

- Display Doors
- Adjustable Shelf

DISPLAY

45682

W : 40, 42, 46, 48
H : 30
D : 14, 17

- Display Doors
- Center Divider
- Adjustable Shelves: 2



**EXPANDED
METAL**



45606

W : 15, 17, 18, 20, 21, 23, 24
H : 30
D : 14, 17

- Right Hinge Framed Expanded Metal Door
- Adjustable Shelf
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**



45607

W : 15, 17, 18, 20, 21, 23, 24
H : 30
D : 14, 17

- Left Hinge Framed Expanded Metal Door
- Adjustable Shelf
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**




45608

W : 30, 34, 36
H : 30
D : 14, 17

- Framed Expanded Metal Doors
- Adjustable Shelf
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**



45609

W : 40, 42, 46, 48
H : 30
D : 14, 17

- Framed Expanded Metal Doors
- Center Divider
- Adjustable Shelves: 2
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

CUBICLE

45720

W : 21, 24
H : 30
D : 14, 17

- 3/4" Adjustable Cubicle Shelves: 4
- Compartments: 6
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE

45722

W : 21, 24
H : 30
D : 14, 17


- Right Hinge Door
- 3/4" Adjustable Cubicle Shelves: 4
- Compartments: 6
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE

45723


W : 21, 24
H : 30
D : 14, 17

- Left Hinge Door
- 3/4" Adjustable Cubicle Shelves: 4
- Compartments: 6
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE

45726


W : 30, 36, 40, 42, 46, 48
H : 30
D : 14, 17

- 3/4" Adjustable Cubicle Shelves: 6
- Compartments: 9
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W), 12 1/4" (40W), 13" (42W), 14 1/4" (46W), 15" (48W)

CUBICLE

45727

W : 30, 36, 40, 42, 46, 48
H : 30
D : 14, 17

- Doors
- 3/4" Adjustable Cubicle Shelves: 6
- Compartments: 9
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W), 12 1/4" (40W), 13" (42W), 14 1/4" (46W), 15" (48W)

CUBICLE

45739

W : 36, 40, 42, 46, 48
H : 30
D : 14, 17

- 3/4" Adjustable Cubicle Shelves: 8
- Compartments: 12
- Compartment Spacing: 3"H min.
- Compartment Widths: 8" (36W), 9" (40W), 9 1/2" (42W), 10 1/2" (46W), 11" (48W)

CUBICLE


45740

W : 36, 40, 42, 46, 48
H : 30
D : 14, 17

- Doors
- 3/4" Adjustable Cubicle Shelves: 8
- Compartments: 12
- Compartment Spacing: 3"H min.
- Compartment Widths: 8" (36W), 9" (40W), 9 1/2" (42W), 10 1/2" (46W), 11" (48W)



**BLIND
CORNER**




45700

W : 33, 36, 39
H : 30
D : 14

W : 36, 39, 42
H : 30
D : 17

- Right Hinge Door
- Adjustable Shelf

**DIAGONAL
CORNER**




45708

W : 25
H : 30
D : 14

W : 28
H : 30
D : 17

- Fixed Shelf

**BLIND
CORNER**



45701

W : 33, 36, 39
H : 30
D : 14

W : 36, 39, 42
H : 30
D : 17

- Left Hinge Door
- Adjustable Shelf

**DIAGONAL
CORNER**




45711

W : 25
H : 30
D : 14

W : 28
H : 30
D : 17

- Right Hinge Door
- Fixed Shelf

**BLIND
CORNER**



45702

W : 33, 36, 39
H : 30
D : 14

W : 36, 39, 42
H : 30
D : 17

- Right Hinge Display Door
- Adjustable Shelf

**DIAGONAL
CORNER**




45712

W : 25
H : 30
D : 14

W : 28
H : 30
D : 17

- Left Hinge Door
- Fixed Shelf

**BLIND
CORNER**



45703

W : 33, 36, 39
H : 30
D : 14

W : 36, 39, 42
H : 30
D : 17

- Left Hinge Display Door
- Adjustable Shelf


DISPENSING


45622

W : 21, 24
H : 30
D : 14, 17

- Right Hinge Door
- Towel Dispenser
- Adjustable Shelf

DISPENSING


45624

W : 18, 21, 24
H : 30
D : 14, 17

- Right Hinge Door
- Towel and Glove/Tissue Dispensers


DISPENSING


45623

W : 21, 24
H : 30
D : 14, 17

- Left Hinge Door
- Towel Dispenser
- Adjustable Shelf


DISPENSING


45625

W : 18, 21, 24
H : 30
D : 14, 17

- Left Hinge Door
- Towel and Glove/Tissue Dispensers


DISPENSING


45626

W : 36
H : 30
D : 14, 17


- Doors
- Towel Dispenser
- Adjustable Shelf

DISPENSING



45627

W : 36
H : 30
D : 14, 17

- Doors
- Towel and Glove/Tissue Dispensers



TOWEL DISPENSER
Bottom-access dispensing of paper towels.



GLOVE/TISSUE DISPENSER
Face-of-door access for gloves or tissues.



SHELF



45601


W : 12, 14, 15, 17, 18, 20, 21, 23, 24, 28,
30, 34, 36

H : 36

D : 14, 17

- Open Storage
- Adjustable Shelves: 2

SHELF



45638

W : 24, 28, 30, 34, 36, 40, 42, 46, 48

H : 36

D : 14, 17

- Center Divider
- Adjustable Shelves: 4

SHELF



45620


W : 12, 14, 15, 17, 18, 20, 21, 23, 24

H : 36

D : 14, 17

- Right Hinge Door
- Adjustable Shelves: 2

SHELF



45681

W : 40, 42, 46, 48

H : 36

D : 14, 17

- Doors
- Center Divider
- Adjustable Shelves: 4

SHELF



45621


W : 12, 14, 15, 17, 18, 20, 21, 23, 24

H : 36

D : 14, 17

- Left Hinge Door
- Adjustable Shelves: 2

SHELF



45650

W : 36, 40, 42, 46, 48

H : 36

D : 14, 17

- Sliding Doors
- Center Divider
- Adjustable Shelves: 4

SHELF



45629

W : 24, 28, 30, 34, 36

H : 36

D : 14, 17

- Doors
- Adjustable Shelves: 2

**MICROWAVE
SHELF**



45649

W : 30

H : 36

D : 14, 17

- Upper Doors
- Microwave Compartment: 16 1/16"H
- 20"D Attachable Microwave Shelf
- Accommodates Typical 1 cu.ft.-2 cu.ft. Microwave

DISPLAY



45630

W : 15, 17, 18, 20, 21, 23, 24
H : 36
D : 14, 17

- Right Hinge Display Door
- Adjustable Shelves: 2

DISPLAY



45654

W : 36, 40, 42, 46, 48
H : 36
D : 14, 17

- Sliding Display Doors
- Center Divider
- Adjustable Shelves: 4

DISPLAY




45631

W : 15, 17, 18, 20, 21, 23, 24
H : 36
D : 14, 17

- Left Hinge Display Door
- Adjustable Shelves: 2

DISPLAY



45658

W : 36, 40, 42, 46, 48
H : 36
D : 14, 17

- Sliding Glass Doors
- Center Divider
- Adjustable Shelves: 4

DISPLAY



45636

W : 30, 34, 36
H : 36
D : 14, 17

- Display Doors
- Adjustable Shelves: 2

DISPLAY



45682

W : 40, 42, 46, 48
H : 36
D : 14, 17

- Display Doors
- Center Divider
- Adjustable Shelves: 4



**EXPANDED
METAL**



45606

W : 15, 17, 18, 20, 21, 23, 24
H : 36
D : 14, 17

- Right Hinge Framed Expanded Metal Door
- Adjustable Shelves: 2
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**



45607

W : 15, 17, 18, 20, 21, 23, 24
H : 36
D : 14, 17

- Left Hinge Framed Expanded Metal Door
- Adjustable Shelves: 2
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**




45608

W : 30, 34, 36
H : 36
D : 14, 17

- Framed Expanded Metal Doors
- Adjustable Shelves: 2
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**



45609

W : 40, 42, 46, 48
H : 36
D : 14, 17

- Framed Expanded Metal Doors
- Center Divider
- Adjustable Shelves: 4
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

CUBICLE

45720

W : 21, 24
H : 36
D : 14, 17

- 3/4" Adjustable Cubicle Shelves: 6
- Compartments: 8
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE

45722

W : 21, 24
H : 36
D : 14, 17


- Right Hinge Door
- 3/4" Adjustable Cubicle Shelves: 6
- Compartments: 8
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE

45723


W : 21, 24
H : 36
D : 14, 17

- Left Hinge Door
- 3/4" Adjustable Cubicle Shelves: 6
- Compartments: 8
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE

45726


W : 30, 36, 40, 42, 46, 48
H : 36
D : 14, 17

- 3/4" Adjustable Cubicle Shelves: 9
- Compartments: 12
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W), 12 1/4" (40W), 13" (42W), 14 1/4" (46W), 15" (48W)

CUBICLE

45727


W : 30, 36, 40, 42, 46, 48
H : 36
D : 14, 17

- Doors
- 3/4" Adjustable Cubicle Shelves: 9
- Compartments: 12
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W), 12 1/4" (40W), 13" (42W), 14 1/4" (46W), 15" (48W)

CUBICLE

45739

W : 36, 40, 42, 46, 48
H : 36
D : 14, 17

- 3/4" Adjustable Cubicle Shelves: 12
- Compartments: 16
- Compartment Spacing: 3"H min.
- Compartment Widths: 8" (36W), 9" (40W), 9 1/2" (42W), 10 1/2" (46W), 11" (48W)


CUBICLE

45740

W : 36, 40, 42, 46, 48
H : 36
D : 14, 17

- Doors
- 3/4" Adjustable Cubicle Shelves: 12
- Compartments: 16
- Compartment Spacing: 3"H min.
- Compartment Widths: 8" (36W), 9" (40W), 9 1/2" (42W), 10 1/2" (46W), 11" (48W)



**BLIND
CORNER**



45700

W : 33, 36, 39
H : 36
D : 14

W : 36, 39, 42
H : 36
D : 17

- Right Hinge Door
- Adjustable Shelves: 2

**DIAGONAL
CORNER**




45708

W : 25
H : 36
D : 14

W : 28
H : 36
D : 17

- Fixed Shelves: 2

**BLIND
CORNER**



45701

W : 33, 36, 39
H : 36
D : 14

W : 36, 39, 42
H : 36
D : 17

- Left Hinge Door
- Adjustable Shelves: 2

**DIAGONAL
CORNER**



45711

W : 25
H : 36
D : 14

W : 28
H : 36
D : 17

- Right Hinge Door
- Fixed Shelves: 2

**BLIND
CORNER**



45702

W : 33, 36, 39
H : 36
D : 14

W : 36, 39, 42
H : 36
D : 17

- Right Hinge Display Door
- Adjustable Shelves: 2

**DIAGONAL
CORNER**



45712

W : 25
H : 36
D : 14

W : 28
H : 36
D : 17

- Left Hinge Door
- Fixed Shelves: 2

**BLIND
CORNER**



45703

W : 33, 36, 39
H : 36
D : 14

W : 36, 39, 42
H : 36
D : 17

- Left Hinge Display Door
- Adjustable Shelves: 2

DISPENSING



45622

W : 21, 24
H : 36
D : 14, 17

- Right Hinge Door
- Towel Dispenser
- Adjustable Shelves: 2

DISPENSING



45624

W : 18, 21, 24
H : 36
D : 14, 17

- Right Hinge Door
- Towel and Glove/Tissue Dispensers
- Adjustable Shelf

DISPENSING



45623

W : 21, 24
H : 36
D : 14, 17

- Left Hinge Door
- Towel Dispenser
- Adjustable Shelves: 2

DISPENSING




45625

W : 18, 21, 24
H : 36
D : 14, 17

- Left Hinge Door
- Towel and Glove/Tissue Dispensers
- Adjustable Shelf

DISPENSING




45626

W : 36
H : 36
D : 14, 17

- Doors
- Towel Dispenser
- Adjustable Shelves: 2


DISPENSING




45627

W : 36
H : 36
D : 14, 17

- Doors
- Towel and Glove/Tissue Dispensers
- Adjustable Shelf



TOWEL DISPENSER
Bottom-access dispensing of paper towels.



GLOVE/TISSUE DISPENSER
Face-of-door access for gloves or tissues.



SHELF



45601

W : 12, 14, 15, 17, 18, 20, 21, 23, 24, 28,
30, 34, 36

H : 42

D : 14, 17

- Open Storage
- Adjustable Shelves: 3

SHELF



45638

W : 24, 28, 30, 34, 36, 40, 42, 46, 48

H : 42

D : 14, 17

- Center Divider
- Adjustable Shelves: 6

SHELF



45620

W : 12, 14, 15, 17, 18, 20, 21, 23, 24

H : 42

D : 14, 17

- Right Hinge Door
- Adjustable Shelves: 3

SHELF



45681

W : 40, 42, 46, 48

H : 42

D : 14, 17

- Doors
- Center Divider
- Adjustable Shelves: 6

SHELF



45621


W : 12, 14, 15, 17, 18, 20, 21, 23, 24

H : 42

D : 14, 17

- Left Hinge Door
- Adjustable Shelves: 3

SHELF



45650

W : 36, 40, 42, 46, 48

H : 42

D : 14, 17

- Sliding Doors
- Center Divider
- Adjustable Shelves: 6

SHELF



45629

W : 24, 28, 30, 34, 36

H : 42

D : 14, 17

- Doors
- Adjustable Shelves: 3

DISPLAY




45630

W : 15, 17, 18, 20, 21, 23, 24
H : 42
D : 14, 17

- Right Hinge Display Door
- Adjustable Shelves: 3

DISPLAY



45654

W : 36, 40, 42, 46, 48
H : 42
D : 14, 17

- Sliding Display Doors
- Center Divider
- Adjustable Shelves: 6

DISPLAY




45631

W : 15, 17, 18, 20, 21, 23, 24
H : 42
D : 14, 17

- Left Hinge Display Door
- Adjustable Shelves: 3

DISPLAY



45658

W : 36, 40, 42, 46, 48
H : 42
D : 14, 17

- Sliding Glass Doors
- Center Divider
- Adjustable Shelves: 6

DISPLAY




45636

W : 30, 34, 36
H : 42
D : 14, 17

- Display Doors
- Adjustable Shelves: 3

DISPLAY




45682

W : 40, 42, 46, 48
H : 42
D : 14, 17

- Display Doors
- Center Divider
- Adjustable Shelves: 6



**EXPANDED
METAL**




45606

W : 15, 17, 18, 20, 21, 23, 24
H : 42
D : 14, 17

- Right Hinge Framed Expanded Metal Door
- Adjustable Shelves: 3
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**




45607

W : 15, 17, 18, 20, 21, 23, 24
H : 42
D : 14, 17

- Left Hinge Framed Expanded Metal Door
- Adjustable Shelves: 3
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**




45608

W : 30, 34, 36
H : 42
D : 14, 17

- Framed Expanded Metal Doors
- Adjustable Shelves: 3
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**



45609

W : 40, 42, 46, 48
H : 42
D : 14, 17

- Framed Expanded Metal Doors
- Center Divider
- Adjustable Shelves: 6
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

CUBICLE

45720

W : 21, 24
H : 42
D : 14, 17

- 3/4" Adjustable Cubicle Shelves: 8
- Compartments: 10
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE

45722

W : 21, 24
H : 42
D : 14, 17


- Right Hinge Door
- 3/4" Adjustable Cubicle Shelves: 8
- Compartments: 10
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE

45723

W : 21, 24
H : 42
D : 14, 17

- Left Hinge Door
- 3/4" Adjustable Cubicle Shelves: 8
- Compartments: 10
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE

45726

W : 30, 36, 40, 42, 46, 48
H : 42
D : 14, 17


- 3/4" Adjustable Cubicle Shelves: 12
- Compartments: 15
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W), 12 1/4" (40W), 13" (42W), 14 1/4" (46W), 15" (48W)

CUBICLE

45727

W : 30, 36, 40, 42, 46, 48
H : 42
D : 14, 17

- Doors
- 3/4" Adjustable Cubicle Shelves: 12
- Compartments: 15
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W), 12 1/4" (40W), 13" (42W), 14 1/4" (46W), 15" (48W)

CUBICLE

45739

W : 36, 40, 42, 46, 48
H : 42
D : 14, 17

- 3/4" Adjustable Cubicle Shelves: 16
- Compartments: 20
- Compartment Spacing: 3"H min.
- Compartment Widths: 8" (36W), 9" (40W), 9 1/2" (42W), 10 1/2" (46W), 11" (48W)

CUBICLE


45740

W : 36, 40, 42, 46, 48
H : 42
D : 14, 17

- Doors
- 3/4" Adjustable Cubicle Shelves: 16
- Compartments: 20
- Compartment Spacing: 3"H min.
- Compartment Widths: 8" (36W), 9" (40W), 9 1/2" (42W), 10 1/2" (46W), 11" (48W)



**BLIND
CORNER**



45700

W : 33, 36, 39
H : 42
D : 14

W : 36, 39, 42
H : 42
D : 17

- Right Hinge Door
- Adjustable Shelves: 3

**DIAGONAL
CORNER**




45708

W : 25
H : 42
D : 14

W : 28
H : 42
D : 17

- Fixed Shelves: 2

**BLIND
CORNER**



45701

W : 33, 36, 39
H : 42
D : 14

W : 36, 39, 42
H : 42
D : 17

- Left Hinge Door
- Adjustable Shelves: 3

**DIAGONAL
CORNER**



45711

W : 25
H : 42
D : 14

W : 28
H : 42
D : 17

- Right Hinge Door
- Fixed Shelves: 2

**BLIND
CORNER**



45702

W : 33, 36, 39
H : 42
D : 14

W : 36, 39, 42
H : 42
D : 17

- Right Hinge Display Door
- Adjustable Shelves: 3

**DIAGONAL
CORNER**




45712

W : 25
H : 42
D : 14

W : 28
H : 42
D : 17

- Left Hinge Door
- Fixed Shelves: 2

**BLIND
CORNER**



45703

W : 33, 36, 39
H : 42
D : 14

W : 36, 39, 42
H : 42
D : 17

- Left Hinge Display Door
- Adjustable Shelves: 3

SHELF

45601

W : 12, 14, 15, 17, 18, 20, 21, 23, 24, 28,
30, 34, 36

H : 47

D : 14, 17

- Open Storage
- Adjustable Shelves: 3

SHELF

45638

W : 24, 28, 30, 34, 36, 40, 42, 46, 48

H : 47

D : 14, 17

- Center Divider
- Adjustable Shelves: 6

SHELF

45620

W : 12, 14, 15, 17, 18, 20, 21, 23, 24

H : 47

D : 14, 17

- Right Hinge Door
- Adjustable Shelves: 3

SHELF

45681

W : 40, 42, 46, 48

H : 47

D : 14, 17

- Doors
- Center Divider
- Adjustable Shelves: 6

SHELF

45621

W : 12, 14, 15, 17, 18, 20, 21, 23, 24

H : 47

D : 14, 17

- Left Hinge Door
- Adjustable Shelves: 3

SHELF


45650

W : 36, 40, 42, 46, 48

H : 47

D : 14, 17

- Sliding Doors
- Center Divider
- Adjustable Shelves: 6

SHELF

45629

W : 24, 28, 30, 34, 36

H : 47

D : 14, 17

- Doors
- Adjustable Shelves: 3



DISPLAY



45630

W : 15, 17, 18, 20, 21, 23, 24
H : 47
D : 14, 17

- Right Hinge Display Door
- Adjustable Shelves: 3

DISPLAY



45654

W : 36, 40, 42, 46, 48
H : 47
D : 14, 17

- Sliding Display Doors
- Center Divider
- Adjustable Shelves: 6

DISPLAY




45631

W : 15, 17, 18, 20, 21, 23, 24
H : 47
D : 14, 17

- Left Hinge Display Door
- Adjustable Shelves: 3

DISPLAY



45658

W : 36, 40, 42, 46, 48
H : 47
D : 14, 17

- Sliding Glass Doors
- Center Divider
- Adjustable Shelves: 6

DISPLAY




45636

W : 30, 34, 36
H : 47
D : 14, 17

- Display Doors
- Adjustable Shelves: 3

DISPLAY



45682

W : 40, 42, 46, 48
H : 47
D : 14, 17


- Display Doors
- Center Divider
- Adjustable Shelves: 6

**EXPANDED
METAL**

45606

W : 15, 17, 18, 20, 21, 23, 24
H : 47
D : 14, 17

- Right Hinge Framed Expanded Metal Door
- Adjustable Shelves: 3
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**

45607

W : 15, 17, 18, 20, 21, 23, 24
H : 47
D : 14, 17


- Left Hinge Framed Expanded Metal Door
- Adjustable Shelves: 3
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**

45608

W : 30, 34, 36
H : 47
D : 14, 17

- Framed Expanded Metal Doors
- Adjustable Shelves: 3
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)

**EXPANDED
METAL**

45609

W : 40, 42, 46, 48
H : 47
D : 14, 17

- Framed Expanded Metal Doors
- Center Divider
- Adjustable Shelves: 6
- Expanded Metal Available in Black, White, Grey (color linked to accent color option)



CUBICLE



45720

W : 21, 24
H : 47
D : 14, 17

- 3/4" Adjustable Cubicle Shelves: 8
- Compartments: 10
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE



45722

W : 21, 24
H : 47
D : 14, 17

- Right Hinge Door
- 3/4" Adjustable Cubicle Shelves: 8
- Compartments: 10
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE




45723

W : 21, 24
H : 47
D : 14, 17

- Left Hinge Door
- 3/4" Adjustable Cubicle Shelves: 8
- Compartments: 10
- Compartment Spacing: 3"H min.
- Compartment Widths: 9 3/8" (21W), 10 7/8" (24W)

CUBICLE



45726

W : 30, 36, 40, 42, 46, 48
H : 47
D : 14, 17

- 3/4" Adjustable Cubicle Shelves: 12
- Compartments: 15
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W), 12 1/4" (40W), 13" (42W), 14 1/4" (46W), 15" (48W)

CUBICLE




45727

W : 30, 36, 40, 42, 46, 48
H : 47
D : 14, 17

- Doors
- 3/4" Adjustable Cubicle Shelves: 12
- Compartments: 15
- Compartment Spacing: 3"H min.
- Compartment Widths: 9" (30W), 11" (36W), 12 1/4" (40W), 13" (42W), 14 1/4" (46W), 15" (48W)

CUBICLE



45739

W : 36, 40, 42, 46, 48
H : 47
D : 14, 17

- 3/4" Adjustable Cubicle Shelves: 16
- Compartments: 20
- Compartment Spacing: 3"H min.
- Compartment Widths: 8" (36W), 9" (40W), 9 1/2" (42W), 10 1/2" (46W), 11" (48W)


CUBICLE



45740

W : 36, 40, 42, 46, 48
H : 47
D : 14, 17

- Doors
- 3/4" Adjustable Cubicle Shelves: 16
- Compartments: 20
- Compartment Spacing: 3"H min.
- Compartment Widths: 8" (36W), 9" (40W), 9 1/2" (42W), 10 1/2" (46W), 11" (48W)

**BLIND
CORNER**

45700

W : 33, 36, 39
H : 47
D : 14

W : 36, 39, 42
H : 47
D : 17

- Right Hinge Door
- Adjustable Shelves: 3


**DIAGONAL
CORNER**

45708

W : 25
H : 47
D : 14

W : 28
H : 47
D : 17

- Fixed Shelves: 2

**BLIND
CORNER**

45701

W : 33, 36, 39
H : 47
D : 14

W : 36, 39, 42
H : 47
D : 17

- Left Hinge Door
- Adjustable Shelves: 3


**DIAGONAL
CORNER**

45711

W : 25
H : 47
D : 14

W : 28
H : 47
D : 17

- Right Hinge Door
- Fixed Shelves: 2

**BLIND
CORNER**

45702

W : 33, 36, 39
H : 47
D : 14

W : 36, 39, 42
H : 47
D : 17

- Right Hinge Display Door
- Adjustable Shelves: 3


**DIAGONAL
CORNER**

45712

W : 25
H : 47
D : 14

W : 28
H : 47
D : 17

- Left Hinge Door
- Fixed Shelves: 2

**BLIND
CORNER**

45703

W : 33, 36, 39
H : 47
D : 14

W : 36, 39, 42
H : 47
D : 17

- Left Hinge Display Door
- Adjustable Shelves: 3

SOFFIT TRIM




15801

W : 96
H : 00
D : 03

- Field Trim to Fit
- Actual Thickness: 3/4"
- Actual Depth: 3"
- Edgebanding: One Long Edge in Casebody Trim Color and 0.5mm Thickness

END SOFFIT CLOSURE



Shown with
15822/15823
Soffit Fascia

15820

W : 00
H : 12, 18, 24
D : 14, 17

- Vertical Grain
- Field Trim to Fit
- Edgebanding: All Edges in Casebody Trim Color and 0.5mm Thickness

VERTICAL SCRIBE




15802

W : 00
H : 96
D : 14, 17

- Field Trim to Fit
- Actual Thickness: 3/4"
- Actual Depth: 3"
- Edgebanding: One Long Edge in Casebody Trim Color and 0.5mm Thickness

SOFFIT FASCIA




Shown with
15820 End
Soffit Closure

15822

W : 12, 14, 15, 17, 18, 20, 21, 23,
24, 28, 30, 34, 36, 40, 42,
46, 48, 96
H : 12, 18, 24
D : 00

- Horizontal Grain
- Field Trim to Fit
- Edgebanding: One Long and Two Short Edges in Casebody Trim Color and 0.5mm Thickness

FRONT FILLER




15803

W : 02, 03, 06, 09, 12
H : 18, 25, 30, 36, 42, 47
D : 14, 17

- Field Trim to Fit
- Includes Top and Bottom Closure Panels
- Edgebanding: Two Short Edges in Casebody Trim Color and 0.5mm Thickness

SOFFIT FASCIA




Shown with
15820 End
Soffit Closure

15823

W : 12, 14, 15, 17, 18, 20, 21, 23,
24, 28, 30, 34, 36, 40, 42,
46, 48
H : 12, 18, 24
D : 00

- Vertical Grain
- Field Trim to Fit
- Edgebanding: One Short and Two Long Edges in Casebody Trim Color and 0.5mm Thickness

CORNER FILLER/CLOSURE




15805

W : 02, 03, 06, 09, 12
H : 18, 25, 30, 36, 42, 47
D : 14, 17

- Field Trim to Fit
- Edgebanding (Corner Filler): Two Short Edges in Casebody Trim Color and 0.5mm Thickness
- Edgebanding (Closure): None

SOFFIT EXTENSION



15824

W : 12, 14, 15, 17, 18, 20, 21, 23,
24, 28, 30, 34, 36, 40, 42,
46, 48
H : 6"-34" (1" increments)
D : 14, 17

- Edgebanding (Box): All Edges in Casebody Trim Color and Thickness
- Edgebanding (Soffit Panel): All Edges in Door/Drawer Trim Color and Thickness

VALANCE
RAIL




15812

W : 48, 96
H : 03
D : 00


- Field Trim to Fit
- Edgebanding: One Long Edge in Casebody Trim Color and Thickness

VALANCE

Advantage Scientific valance and slope valance units are designed to allow continuous valance aprons under wall units. A finished end/drop left or drop right indication lowers the designated end panel to align with the valance rail. Field cutting of intermediate sides is not required to install under-cabinet valance lighting. Return corners are not required upon installation.



LEFT FRONT CORNER



Because cabinet bottom is set up in relation to door bottom edge, exposed hinges will not align with conventional cabinet hinge locations. Pulls and locks will not align with conventional pull and lock locations.