

APPARATUS COMPONENTS

2

RODS 2
CLAMPS 2
SOCKETS 2

REV. 06/13/24

©2023 STEVENS INDUSTRIES, INC.
ALL RIGHTS RESERVED.
Due to continuing product
development, Stevens reserves
the right to change specifications
without notice.

3/4" HORIZONTAL APPARATUS RODS



522H1-24

WIDTH

24"	522H1-24
30"	522H1-30
36"	522H1-36
48"	522H1-48

Satin aluminum finish
3/4" horizontal rods with
rounded ends.

3/4" VERTICAL APPARATUS RODS



522V1-24

HEIGHT

24"	522V1-24
30"	522V1-30
36"	522V1-36

Satin aluminum finish
3/4" vertical rods with
top end rounded and
bottom end tapered to
fit into burette socket.

3/4" ROD CLAMP



522C1

Satin aluminum
finish closed style
clamp for holding
rods perpendicular to
each other. Adjustable
and furnished with
set screws. For 3/4"
diameter horizontal and
vertical rods.

3/4" RECESSED BURETTE SOCKET



For Solid Core Tops

522RB

Satin aluminum finish
rod socket with tapered
bore and flat top for
recess mounting in
countertop. For 3/4"
diameter vertical rods.

3/4" SURFACE MOUNT BURETTE SOCKET



For Laminate Tops

522SB

Satin aluminum finish
rod socket with tapered
bore and radiused top
for surface mounting
on countertop. For 3/4"
diameter vertical rods.