

PEGBOARDS

2

ADVANTAGE SCIENTIFIC
PEGBOARDS 2


ADVANTAGE SCIENTIFIC
STAINLESS STEEL
PEGBOARDS 2

REV. 6/20/23

©2023 STEVENS INDUSTRIES, INC.
ALL RIGHTS RESERVED.
Due to continuing product
development, Stevens reserves
the right to change specifications
without notice.

ADVANTAGE SCIENTIFIC PEGBOARDS

ADVANTAGE SCIENTIFIC PEGBOARDS




- 1" Material
- Standard Color: Black
- Special Order Colors: Charcoal, White
- Black Pegs and Mounting Screws Included
- Stainless Steel Drip Tray Available Separately

STOCKED ITEM	ADVANTAGE SCIENTIFIC PEGBOARDS				STAINLESS STEEL DRIP TRAY (AVAILABLE SEPARATELY)
	MODEL #	W	H	# PEGS	
	52D1815	18	15	15	52T184
	52D1824	18	24	21	52T184
•	52D2415	24	15	20	52T244
	52D2418	24	18	20	52T244
	52D2424	24	24	28	52T244
	52D2430	24	30	36	52T244
•	52D2436	24	36	40	52T244
	52D3024	30	24	40	52T304
	52D3030	30	30	50	52T304
•	52D3036	30	36	55	52T304
	52D3230	32	30	55	52T324
•	52D3236	32	36	61	52T324
•	52D3624	36	24	42	52T364
•	52D3630	36	30	60	52T364
	52D3636	36	36	66	52T364

ADVANTAGE SCIENTIFIC STAINLESS STEEL PEGBOARDS

ADVANTAGE SCIENTIFIC STAINLESS STEEL PEGBOARDS



- 20 Gauge Type 304 Stainless Steel
- Black Pegs and Mounting Screws Included
- Stainless Steel Drip Tray Available Separately

STOCKED ITEM	ADVANTAGE SCIENTIFIC STAINLESS STEEL PEGBOARDS				STAINLESS STEEL DRIP TRAY (AVAILABLE SEPARATELY)
	MODEL #	W	H	# PEGS	
	52S1815	18	15	15	52T184
	52S1824	18	24	21	52T184
•	52S2415	24	15	20	52T244
	52S2418	24	18	20	52T244
	52S2424	24	24	28	52T244
	52S2430	24	30	36	52T244
•	52S2436	24	36	40	52T244
	52S3024	30	24	40	52T304
	52S3030	30	30	50	52T304
•	52S3036	30	36	55	52T304
	52S3230	32	30	55	52T324
•	52S3236	32	36	61	52T324
•	52S3624	36	24	42	52T364
•	52S3630	36	30	60	52T364
	52S3636	36	36	66	52T364

ADVANTAGE SCIENTIFIC STAINLESS STEEL DRIP TRAYS

STAINLESS STEEL DRIP TRAYS



- Plastic Drain Tubing Included

STOCKED ITEM	MODEL #	W	H	D
•	52T184	18 1/4	1	4
•	52T244	24 1/4	1	4
•	52T304	30 1/4	1	4
•	52T324	32 1/4	1	4
•	52T364	36 1/4	1	4

