

EPOXY SINKS FOR EPOXY TOPS	3
EPOXY SINKS FOR PHENOLIC TOPS	6
SOLID SINK ACCESSORIES	8
STAINLESS STEEL SINKS	9
STAINLESS STEEL SINK ACCESSORIES	12

Rev. 08-27-24


©2023 STEVENS INDUSTRIES, INC.
ALL RIGHTS RESERVED.
Due to continuing product
development, Stevens reserves
the right to change specifications
without notice.





**ADVANTAGE
SCIENTIFIC**
BY STEVENS

EPOXY SINKS FOR EPOXY TOPS

DROP-IN RECTANGULAR EPOXY SINKS - BLACK





Model #	Outlet	Dimensions (in)		
		Inside (nominal)		
		A	B	C
52D03	Center	12	8	6
52D05	Center	14	10	6
52D10	Center	16	8	7
52D15	Center	16	12	8
52D20	Corner	16	16	8
52D22	Center	18	7	6
52D25	Corner	18	15	8
52D30	Corner	18	15	11
52D54	Corner	25	15	8
52D55	Corner	25	15	10
52D59	Corner	28	15	12




ADA DROP-IN RECTANGULAR EPOXY SINKS - BLACK




Model #	Outlet	Dimensions (in)		
		Inside (nominal)		
		A	B	C
52D05A	Corner	14	10	5
52D07A	Corner	14	14	5
52D25A	Corner	18	15	5
52D26A	Corner	18	15	5/11
52D55A	Corner	25	15	5
52D56A	Corner	25	15	5/11



DROP-IN TROUGH EPOXY SINKS - BLACK




Model #	Outlet	Dimensions (in)		
		Inside (nominal)		
		A	B	C
52DT36	Off-Center	35	6	4
52DT48	Off-Center	47	6	4





SINKS FOR EPOXY TOPS

DROP-IN HEXAGONAL EPOXY SINK




52DH20

Outlet	Dimensions (in)		
	Inside (nominal)		
	L	W	D
Center	30	26	7




OVAL CUP EPOXY SINK



52CS4R

Dimensions (in)					
Inside			Outside		
L	W	D	L	W	D
5.7	2.7	3.7	7.5	4.5	8.3


END AND WALL EPOXY SINK



52E46

Outlet	Dimensions (in)			
	Inside			Back Splash
	L	W	D	
Corner	22.0	13.6	12.6	3.8

STAND FOR 52E46 SINK



49739

Dimensions (in)		
L	W	D
20.5	16.1	13.75

EPOXY SINK ACCESSORIES FOR EPOXY TOPS

OPEN END
OVERFLOW
- BLACK



52POB

Model #	Dimensions (in)	
	A	B
52POB4	1 5/8	4
52POB6	1 5/8	6
52POB8	1 5/8	8




OUTLET/STOPPER
- BLACK



52P98

Dimensions (in)		
Size	A	B
1 1/2 NPSL	3 3/8	3 3/4





EPOXY SINKS FOR PHENOLIC TOPS

DROP-IN RECTANGULAR EPOXY SINKS - BLACK



Phenolic Sinks are shipped with Outlet, Stopper and Overflow.

Model #	Outlet	Dimensions (in)		
		Inside (nominal)		
		A	B	C
52S03	Corner	12	8	8
52S05	Corner	14	10	6
52S10	Corner	16	8	7
52S15	Corner	16	12	8
52S25	Corner	18	15	8
52S30	Corner	18	15	11
52S55	Corner	25	15	10
52S59	Corner	28	15	12



ADA DROP-IN RECTANGULAR EPOXY SINKS - BLACK




Phenolic Sinks are shipped with Outlet, Stopper and Overflow.

Model #	Outlet	Dimensions (in)		
		Inside (nominal)		
		A	B	C
52S25A	Corner	18	15	5
52S55A	Corner	25	15	5


DROP-IN HEXAGONAL EPOXY SINKS



Model #	Outlet	Dimensions (in)		
		Inside (nominal)		
		L	W	D
52SH20	Offset	30	26	7




DROP-IN TROUGH EPOXY SINKS - BLACK



See Epoxy Sink Accessories Page for Outlet, Stopper and Overflow.

Model #	Outlet	Dimensions (in)		
		Inside (nominal)		
		A	B	C
52DT36	Off-Center	35	6	4
52DT48	Off-Center	47	6	4



EPOXY SINK ACCESSORIES FOR PHENOLIC TOPS

OPEN END OVERFLOW - BLACK




Phenolic Sinks are shipped with Outlet, Stopper and Overflow, except Trough Sinks.

Model #	Dimensions (in)	
	A	B
52SOT4	1 1/2	4
52SOT6	1 1/2	6
52SOT8	1 1/2	8




OUTLET/STOPPER - BLACK




Phenolic Sinks are shipped with Outlet, Stopper and Overflow, except Trough Sinks.

Model #	Dimensions (in)		
	Size	A	B
52SOS1	1 1/2 NPSL	2 7/8	4 5/32



GENERAL SOLID TOP SINK ACCESSORIES

SINK COVER - MATCHES TOP COLOR




Phenolic construction for phenolic and epoxy rectangular and trough drop-in sinks. Bottom extension is sized to fit sink inside dimensions.

Model #	Sink Size (in)	
	Inside (nominal)	
	A	B
52C03	12	8
52C05	14	10
52C10	16	8
52C15	16	12
52C20	16	16
52C22	18	6.5
52C25	18	15
52C30	18	15
52C50	24	16
52C54	25	15
52C55	25	15
52C59	28	15
52C36	35	6
52C48	47	6



SIDE PROFILE

ALL WAY TRAP



52P91

Can be used in "P" or "S" configuration.

Inlet	Outlet	Dimensions (in)						
		A	A1	B		B1		C
				MIN	MAX	MIN	MAX	
1 1/2 NPS	1 1/2 VMT-VS	9 1/2	6	2 1/2	7 1/4	6 1/8	10 7/8	2 1/2

STAINLESS STEEL SINKS

CUSTOM SINGLE BOWL

18-GA./304 ALLOY



Must specify faucet hole punching required. (See page 10)

Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
52SL1515AGR	15	6	15	12	9	3 1/2
52SL17519AGR	19	7 1/2	17 1/2	16	11 1/2	3 1/2
52SL1815AGR	15	7 1/2	18	12	12	3 1/2
52SL1921AGR	21	7 1/2	19	18	14	3 1/2

CUSTOM DEEP SINGLE BOWL

18-GA./304 ALLOY



Must specify faucet hole punching required. (See page 10)

Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
52SLX1815AGR	15	9	18	12	12	3 1/2
52SLX1921AGR	21	10 1/2	19	18	14	3 1/2

CUSTOM ADA SINGLE BOWL

18-GA./304 ALLOY



Must specify faucet hole punching required. (See page 10)

Specify depth in 1/2" increments, from 4.5" to 6.5".

Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H specify depth in 1/2" increments (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
52SLADA17519AGR	19	4 1/2 - 6 1/2	17 1/2	16	11 1/2	3 1/2
52SLADA1815AGR	15	4 1/2 - 6 1/2	18	12	12	3 1/2
52SLADA1921AGR	21	4 1/2 - 6 1/2	19	14	18	3 1/2

STAINLESS STEEL SINKS

CUSTOM DOUBLE BOWL

18-GA./304 ALLOY



Must specify faucet hole punching required. (See page 10)

Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
52DL1625AGR	25	7 1/2	16	10	10	3 1/2
52DL17537AGR	37	7 1/2	17 1/2	16	11 1/2	3 1/2
52DL1829AGR	29	7 1/2	18	12	12	3 1/2
52DL1925AGR	25	7 1/2	19	10	14	3 1/2
52DL1933AGR	33	7 1/2	19	14	14	3 1/2

CUSTOM DEEP DOUBLE BOWL

18-GA./304 ALLOY



Must specify faucet hole punching required. (See page 10)

Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
52DLX1829AGR	29	9	18	12	12	3 1/2
52DLX1933AGR	33	10 1/2	19	14	14	3 1/2

CUSTOM ADA DOUBLE BOWL

18-GA./304 ALLOY



Must specify faucet hole punching required. (See page 10)

Specify depth in 1/2" increments, from 4.5" to 6.5".

Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H specify depth in 1/2" increments (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
52DLADA17537AGR	37	4 1/2 - 6 1/2	17 1/2	16	11 1/2	3 1/2
52DLADA1829AGR	29	4 1/2 - 6 1/2	18	12	12	3 1/2
52DLADA1925AGR	25	4 1/2 - 6 1/2	19	10	14	3 1/2
52DLADA1933AGR	33	4 1/2 - 6 1/2	19	14	14	3 1/2

STAINLESS STEEL SINKS

SINGLE BOWL

22-GA. STAINLESS STEEL



Must specify faucet hole punching required. (See page 10)

Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
5211515	15	5 1/8	15	12	10	2
5211719	17	6 1/8	19	14	14	3 1/2

DOUBLE BOWL

22-GA. STAINLESS STEEL



Must specify faucet hole punching required. (See page 10)

Model #	Overall Sink Size (in)			Bowl Size (in)		Outlet
	W (left to right)	H (bowl depth)	D (front to back)	W (left to right)	D (front to back)	
5222519	25	6 1/4	19	10	14	3 1/2
5223319	33	6 1/4	19	14	14	3 1/2

Faucet Hole Punching Options – Single Bowl Sinks


Specify faucet hole punching required by adding the appropriate number code.

- 1 for (1) hole, centered
- 2 for (2) holes on 4" centers
- 3 for (3) holes on 4" centers


Example:

Specify 52SL1921AGR-3 for 52SL1921AGR sink with 3 holes on 4" centers


(1) hole, centered



(2) holes on 4" centers



(3) holes on 4" centers



If ADA Sink, will also require depth information. Provide after hole notation.

Example:

52SLADA1815AGR-2-4.5 is a 52SLADA1815AGR sink with 2 holes on 4" centers that is 4.5" deep.

Faucet Hole Punching Options – Double Bowl Sinks


Specify faucet hole punching required by adding the appropriate number code.

- 1 for (1) hole, centered
- 2 for (2) holes on 4" centers
- 3 for (3) holes on 4" centers
- 4 for (4) holes on 4" centers


Example:

Specify 52DL1829AGR-3 for 52DL1829AGR sink with 3 holes on 4" centers


(1) hole, centered




(2) holes on 4" centers



(3) holes on 4" centers



(4) holes on 4" centers



STAINLESS STEEL SINK ACCESSORIES


STRAINER



52J35

Stainless steel cup strainer with removable basket. Basket fitted with rubber stopper. For 3 1/2" sink outlet.

UNDERMOUNT SINK SUPPORT ASSEMBLY



10942

Includes (2) Channels. Channel Sized To Fit Inside 48" Wide Unit.

W: 46
H: 03
D: 02


DRAIN WITH FLAT STRAINER



52J35FS

Chrome plated brass drain with flat strainer. For 3 1/2" sink outlet.

PLASTER AND SEDIMENT TRAP



52P97

Compatible with All Craft Sinks
Black
1 1/2-Gal. Capacity
(2) Inlets (1 1/2")
(1) Outlet (2")
(1) Cleanout (3")

W: 09
H: 12
D: 09

Dimensions (in)

A	B	C
9	12	6 7/8


STRAINER



52D5018A

Stainless steel cup strainer with removable basket. For 2" sink outlet.

ALL WAY TRAP



52P91

Can be used in "P" or "S" configuration.

Inlet	Outlet	Dimensions (in)						
		A	A1	B		B1		C
				MIN	MAX	MIN	MAX	
1 1/2 NPS	1 1/2 VMT-VS	9 1/2	6	2 1/2	7 1/4	6 1/8	10 7/8	2 1/2